6.1 FOOD: WHERE DOES IT COME FROM?



How to prepare sprouts from whole seeds?

What is required?

Dry whole seeds of gram/
moong, petridish and
cotton wool.

How will you proceed?

- 1. Wash and soak the seeds in water for a day.
- 2. Drain excess water and leave the soaked seeds overnight in a petridish and cover the seeds with wet cotton.
- 3. Keep the cotton wet (it should not get dry).
- 4. After 2-3 days, do you see any new white structures emerging out of the seeds? Name the process involved?



(a)



(b)



(c)



(d)

Fig.6.1.1Germination of gram seeds

What have you learnt?

The seeds germinate after being kept in moist condition. This is the way a new plant develops from a seed.



Extension

- Do all types of seeds form sprouts?
 Name the conditions required for sprouting.
- 3. Use the sprouts of moong and gram to prepare a nutritious and tasty snack.

Notes			

_ UPPER PRIMARY SCIENCE KIT MANUAL



6.2 Components of Food

ACTIVITY 6.2.1

How can we identify the major components in our food?

How will you proceed?

A. Test for starch

- 1. Take a test tube and put some crushed potatoes in it.
- 2. Add few drops of iodine solution in the test tube.

What is required?

Crushed vaw potatoes,
powdered dal, crushed
groundnuts, dilute
iodine solution, copper
sulphate solution,
sodium hydroxide
solution (caustic soda),
piece of white paper, test
tubes and a dropper.

Observation

You will see dark blue black colour. This indicates the presence of starch in potato.

B. Test for protein

- 1. Take a test tube and put some powdered dal in it.
- 2. Add 2 drops of copper sulphate solution and 10 drops of caustic soda and shake the test tube.

Observation

You will see a violet colour. This indicates the presence of proteins in dal.





C. Test for fats

1. Take crushed ground nuts in between the folds of a white paper and press them with the help of thumb.

Observation

You will see a translucent patch on the paper at the place where the nuts were pressed.

This indicates the presence of fats (oil).

What have you learnt?

The major components of our food are carbohydrates (Starch), proteins and fats.

Extension

Repeat the experiment with powdered rice, powdered gram (besan) or soybean seeds, etc., and fill in the table given below by putting a tick ($\sqrt{\ }$) in the correct place:

Food items	Starch	Protein	fat
Powdered Rice			
Powdered gram			
Powdered soybean seeds			

4 _____UPPER PRIMARY SCIENCE KIT MANUAL





6.3 FIBRE TO FABRIC

ACTIVITY 6.3.1

How do you identify a fibre?

What is required?

Different types of fibres and a magnifying glass.

How will you proceed?

1. Take three types of yarns of equal size, for example a piece of woollen thread, a piece of cotton thread and a piece of synthetic fibre.

2. Spread them on a sheet of white paper. Observe the yarns through the magnifying glass.

3. Identify the yarn on the basis of size and number of 'burrs'. [See Fig. 6.3.1] below.




Fig.6.3.1 Different types of fibres

What have you learnt?

Fibres can be identified by their textures. The 'burrs' in wool are the maximum, less in cotton yarn and the least in synthetic yarn. We can use a magnifying glass to see the burrs clearly.



		•	
Exte	ns	LO	п

1. Identification of different yarns by their smell, and the left-over after burning. (burn the yarns very carefully).

Sample	Smell	Characteristics of left-over after burning
A		
В		
С		

2. Check for absorption of water by different types of fibres.

Notes			

6 _____UPPER PRIMARY SCIENCE KIT MANUAL



6.4 SORTING MATERIALS INTO GROUPS

ACTIVITY 6.4.1

How do you classify materials/objects on the basis of their appearance?

What is required?

Sand paper, card

board, piece of wood,

copper wire,

aluminium sheet,

chalk, medicine

wrappers, iron nails,

zinc pieces and

magnesium ribbon.

How will you proceed?

- 1. Observe each of the above materials and note which of them have shiny surface.
- 2. Take an iron nail and rub it properly with the sand paper. Does the freshly rubbed surface appear shiny?
- 3. Repeat the same procedure with other materials. Why do we need to rub the surface of materials?
- 4. Classify these materials as *lustrous* or *non-lustrous materials* i.e., shiny and non-shiny in appearance.

Sl. No	Object	Lustrous (shiny)	Non–lustrous (non-shiny)
1.			
2.			
3.			
4.			
5.			
6.			
7.			
8.			



- 1. Freshly cut/rubbed surfaces of some materials have shiny appearance i.e. they have lustre.
- 2. The materials which have lustre are copper wire, aluminium sheet, medicine wrappers, iron nails, zinc pieces and magnesium ribbon and those without lustre are card board, piece of wood and chalk.

Extension


Try to classify some more materials/objects on the basis of their appearance.

Notes			

8 ______ UPPER PRIMARY SCIENCE KIT MANUAL







How will you proceed?

- 1. Take some test tubes.
- 2. Fill each of them up to 1/4th with water.
- 3. Add a spatula full of any one substance in one of the test tubes and shake it. Observe what happens?
- 4. Similarly, add small amounts of the other substances in the remaining test tubes and shake them.
- 5. Which substances disappear or dissolve in water and which substances do not disappear or dissolve?
- 6. Separately list the substances which dissolve in water and the ones which do not dissolve.





Fig.6.4.2 Finding solubility of different substances in water

SORTING MATERIALS INTO GROUPS _____



- 1. Some substances dissolve in water while others do not.
- 2. Substances like sugar and salt which disappear or dissolve in water, are called water-soluble.
- 3. Substances like chalk powder, sand and saw dust which do not disappear / dissolve in water are called water-insoluble.


Extension

Notes

Check the solubility of glucose, starch, refined flour, iron filings etc. in water and classify them as water soluble or insoluble.



LUPPER PRIMARY SCIENCE KIT MANUAL



ACTIVITY 6.4.3


What happens when other liquids are added to water?

What is required?

Vinegar, Lemon juice,
alcohol, mustard oil,
coconut oil, cooking oil,
kerosene, test tubes,
dropper, wash bottle.

How will you proceed?

- 1. Take few test tubes.
- 2. Take equal volume of water in each of them i.e., fill each of them approximately 1/3rd with water.
- 3. Add a few drops of oil with a dropper in one of the test tubes [Fig.6.4.3 (a)].
- 4. Similarly, add different liquids to different test tubes containing water [Fig.6.4.3 (b)].
- 5. Shake the test tubes well and observe them after sometime.
- 6. You can see that the liquids either mix with water well or do not mix and form a separate layer [Fig. 6.4.3(a) and (b)].






Fig.6.4.3 (a) liquids which mix well with water

(b) Liquid which do not mix with water

SORTING MATERIALS INTO GROUPS

_ 11



7. Make a list of the liquids which mix well with water and others which do not mix with water.
What have you learnt?

Some liquids completely mix with water while some others do not. Liquids like vinegar, lemon juice and alcohol that mix with water are soluble in water and liquids that do not mix with water like mustard oil, cooking oil, coconut oil and kerosene are said to be insoluble in water.

Notes			

12 _____ UPPER PRIMARY SCIENCE KIT MANUAL



ACTIVITY 6.4.4

Can you see through materials?

What is required?

A piece of wood, a piece of glass, different plastic sheets, a piece of cardboard, white paper, oil, cotton, pencil, a

flower, source of light,

bulb, iron nails.

test tube, light or electric

How will you proceed?

- 1. Take each of the above materials one by one and try to see through them.
- 2. You can do this by placing a pencil or an iron nail behind each one of these materials and try to see it.
- 3. Can you see through them?
- 4. List the materials through which you could see and through which you could not see.
- 5. Take a white paper.

6. Hold the iron nail behind the paper and try to see through it.



Fig.6.4.4 Distinguishing among opaque, transparent and translucent object

SORTING MATERIALS INTO GROUPS _____



- 7. Are you able to see the iron nail through the paper?
- 8. Now put 2-3 drops of oil and smear with the help of a cotton plug so that the oil occupies a larger area.
- 9. Now try to see the nail through this paper.
- 10. What do you see? Can you now see the nail?

- 1. Some materials through which we cannot see are called **opaque** while other materials through which we can see clearly are called **transparent**.
- 2. Some other materials through which one cannot see clearly, are called **translucent**.
- 3. Oiled paper is translucent and we can see partially through the translucent materials.

Extension

Take some more materials like coloured glass, plastic sheets, coloured paper sheets etc., and classify them as opaque, transparent or translucent.


14 ______ UPPER PRIMARY SCIENCE KIT MANUAL







6.5 SEPARATION OF SUBSTANCES



How will you proceed?

1. Take some muddy water in a beaker.

2. Leave it undisturbed for some time till the mud settles down [Fig.6.5.1 (a)]. What

is this process called?

3. Pour the water above mud by tilting the beaker gently into another beaker [Fig.6.5.1 (b)]. What is this process called? Is the water poured in second beaker clear?






Fig.6.5.1 (a) Sedimentation of mud (b) Decantation of mud (c) Filterationh of mud

- If no, what will you do to get clear water? You can proceed as given below.
- 4. Take a fine cloth and gently pour the water through it into a third beaker [Fig. 6.5.1 (c)]

 Is the water which is obtained in the third beaker clear?
- 5. If not, what will you suggest to get the clear water? You can proceed further as given in step 6.
- 6. Pour the water through a filter paper folded in the form of a cone kept in a funnel [Fig. 6.5.1 (d)].





Fig.6.5.1 (d) Filtration of mud from muddy water

How to fold the filter paper? You can follow the steps which are shown below [Fig. 6.5.1 (e)].



UPPER PRIMARY SCIENCE KIT MANUAL



In a mixture of water and a solid, the heavier particles get settled. The water above the settled particles can be removed. In order to further remove the finer particles, a cloth or a filter paper can be used. Mud can be separated from muddy water by sedimentation, decantation and filtration to give the clear water.

Extension

Notes

Think of some more mixtures where these techniques can be used for separation.

1000			

SEPARATION OF SUBSTANCES





How will you proceed?

1. Take a spatula full of common salt and transfer it into a beaker.




Fig.6.5.2 Separation of salty from salty water

- 2. Dissolve it in the minimum amount of water to prepare the salt solution.
- 3. Take a small quantity of this solution in a China dish.

8 _____UPPER PRIMARY SCIENCE KIT MANUAL

- 4. Heat the China dish containing salt solution on a kerosene burner [Fig.6.5.2(a)].
- 5. Do you see any change in the volume of solution?
- 6. Keep heating till all the water evaporates.
- 7. What is left behind in the China dish? [Fig.6.5.2(b)]. What about its quantity?

- 1. When the salt solution is heated, water gets evaporated leaving behind the salt.
- 2. Salt from the salt solution can be obtained by evaporation of water.

Extension


Think of some more substances which can be obtained by evaporation of their solutions. Repeat this activity with these substances also.

SEPARATION OF SUBSTANCES

19










How will you proceed?

- 1. Take a mixture of sand and salt in a beaker and add some water. Why is water added to the mixture?
- 2. Now stir the contents with a glass rod [Fig.6.5.3 (a)].
- 3. Leave them undisturbed for some time. What do you see? [Fig. 6.5.3 (b)].
- 4. Decant the liquid from the beaker to a China dish [Fig.6.5.3 (c)].
- 5. Heat the liquid in the China dish over a kerosene burner. The water will start evaporating [Fig. 6.5.3 (d)].
- 6. What will be left in the China dish after all the water has evaporated?



_____UPPER PRIMARY SCIENCE KIT MANUAL






Fig.6.5.3 (a) Mixture of sand water


- (b) Sedimentation of sand
- (c) Decantation of salty water
- (d) Evaporation of water
- (e) Salt after evaporation

- 1. Salt is soluble in water whereas sand is insoluble.
- 2. When water is added to the mixture of sand and salt, sand being heavier then water, settles down.
- 3. The mixture of sand and salt can be separated by decantation and evaporation.

SEPARATION OF SUBSTANCES

_ 2′





What is required?

Salt, water, beaker,

spatula, glass rod.

How will you proceed?

- 1. Take a 50 mL beaker and half fill with water.
- 2. Add a spatula full of salt to it and stir it well.
- 3. Does the salt disappear?
- 4. Add little more of the salt into it and stir again. What happens now?
- 5. Keep on adding the salt with continuous stirring till the salt stops disappearing (dissolving).
- 6. What will happen if you add some more salt? Does it dissolve?




Fig.6.5.4 Salt dissolves in water at room temperature

22 _____UPPER PRIMARY SCIENCE KIT MANUAL





A given amount of water dissolves only a definite amount of a substance (salt) at a given temperature. The solution so obtained is called a saturated solution.

Extension


- 1. Add some water to the saturated solution containing some undissolved salt and stir.
- 2. What do you observe?

Notes			

SEPARATION OF SUBSTANCES

23





How can a saturated solution be made to dissolve some more substance?

What is required?
Salt, beaker,
spatula, glass rod,
kerosene burner.

How will you proceed?

- 1. In a beaker take the saturated solution of salt as prepared in the previous activity.
- 2. Is there some way that saturated solution could be made to dissolve more salt?
- 3. Add a small quantity of salt and heat the beaker over the kerosene burner [Fig.6.5.5]
- 4. What do you observe? Has the added salt disappeared?
- 5. If you see that the added salt has disappeared, add some more salt.
- 6. What happens now?




Fig.6.5.5 Amount of dissolvable salt increases in hot water

24 ______ UPPER PRIMARY SCIENCE KIT MANUAL



Some more amount of a substance can be made to dissolve in the saturated solution by heating it.

Extension

Try to dissolve the copper sulphate into its saturated solution.

Notes			

SEPARATION OF SUBSTANCES









Are different substances dissolved in varying quantities in water?

What is required?

Common salt, sugar,

spatula, test tubes.

How will you proceed?

- 1. Take a test tube and fill it up to $1/4^{th}$ with water.
- 2. With the help of a spatula, add common salt to the test tube.
- 3. Shake the test tube well to dissolve the salt.
- 4. Keep adding salt using a spatula till you get a saturated solution. Also count the number of times the salt is added to get the saturated salt solution.
- 5. Repeat the above procedure using sugar.
- 6. Are the amounts of salt and sugar required to form their saturated solutions the same?

What have you learnt?

Different substances dissolve to different extent in equal volumes of water.





Fig.6.5.6 Different substances dissolved in varying quantities in equal volume of water

_____ UPPER PRIMARY SCIENCE KIT MANUAL





6.6 CHANGES AROUND US



How will you proceed?

Part A

- 1. Take a rubber band. Observe its shape and size [Fig.6.6.1(a)].
- 2. Can you change its shape or size?
- 3. Stretch it. What do you observe? [Fig.6.6.1(b)].
- 4. Now release it. Does it regain its shape and size?




Fig.6.6.1 (i) Reversable change

Part B

- 1. Now take a balloon. Observe its shape and size.
- 2. Can you change its shape and size? How?
- 3. Inflate the balloon. Observe its shape and size.
- 4. Now release the air from the balloon. Does it regain its original shape and size?






Fig.6.6.1 (ii) Reversible change

(a)

We can get back the same shape and size of some objects. So we understand that some changes can be reversed.

Extension

- 1. Observe the balloon–seller making different types of toys, e.g., cat, monkey etc., from the big balloon.
- 2. Do you observe any change in size only or shape only or both of these in the balloon from which the toy is made?

28 _____ UPPER PRIMARY SCIENCE KIT MANUAL





ACTIVITY 6.6.2

Can you explain the difference between changes which can be reversed and those which cannot be reversed?

How will you proceed?

of same size, scissors. 1. Take a sheet of paper and fold it to make an aeroplane. Do you observe any change in its shape and size?

2. Now unfold it and observe its shape and size. Does the paper regain its original shape and size?

What is required?

Two square sheets of paper





Fig.6.6.2 (i) A toy aeroplane made by folding paper

- 3. Take the other sheet of paper and cut it into four equal parts.
- 4. Observe the size of all the parts.
- 5. Can you get back the original size of sheet of paper again?



CHANGES AROUND US _








Fig.6.6.2 (ii) (a) Plane sheet of paper (b) Sheet divided into four equal parts

Some changes can be reversed while some others cannot be reversed.

Extension

Try to make different paper toys by folding and cutting and then observe which type of change is reversed.

30 ______UPPER PRIMARY SCIENCE KIT MANUAL





Are 'burning of candle and 'melting of wax'

What is required? A candle, a match box, a test tube, some wax and a scale.

How will you proceed?

Part A

- 1. Take a candle and measure its length with a scale Fig.6.6.3 (i).
- 2. Now light the candle and let it burn for sometime.
- 3. Can this change be reversed?
- 4. Now measure the length of the candle again. Can you get the candle of same length again?



Fig.6.6.3 (i) Non-reversible change

Part B

- 1. Take some wax in the test tube.
- Heat the wax in the test tube till it melts completely [Fig.6.6.3 (ii)].
- 3. Now let the wax cool down.
- 5. What do you observe?

CHANGES AROUND US ___






Fig.6.6.3 (ii) Reversible change

Notes

The candle of same length cannot be obtained again i.e., change in length of candle cannot be reversed.

Hence, burning of candle is not reversible but melting of wax is reversible.

10003					

32 _____UPPER PRIMARY SCIENCE KIT MANUAL



6.7 GETTING TO KNOW PLANTS

ACTIVITY 6.7.1

Why does the colour of white petals of a flower change when a flowering twig is dipped in some coloured solution?

How will you proceed?

- Fill one third of a tumbler with coloured (red/blue) water.
- 2. Give an oblique cut to a white flowering twig and put it in the beaker for 3 4 hours.
- 3. Keep this set-up in an open space in bright sunlight.
- 4. Do you observe any change in the colour of the petals of the flower?

What is required?

A white flowering twig of balsam or chandini plant, coloured ink (blue/red), a tumbler, a magnifying glass and

a blade.




Fig.6.7.1 (i) Water moves upwards through the narrow tubes inside the stem of the plants



- 5. Is there any resemblance between the colour of the petals and the solution?
- 6. Now cut the twig across and observe the colour inside with the help of a magnifying glass.

The stem of the twig contains narrow tubes (special cells) for transportation of water and minerals to other parts of the plant.

Extension

Take another white flowering twig and split it halfway along its length and put one of the split ends in tumbler containing red colured solution and the other end into the tumbler containing blue colured solution, containing red (a) and blue (b) coloured solutions.

What do you observe after 3-4 hours?




Fig.6.7.1 (ii) Water moves upwards through special tissues (xylem) of the plants

UPPER PRIMARY SCIENCE KIT MANUAL

34



ACTIVITY 6.7.2

How to Study the venation patterns in leaves?

What is required?

Leaves of plants like
peepal, rose, banana,
bamboo, kaner, wheat,
maize, mango, lemon,
balsam, etc., a magnifying
glass and old newspapers.

How will you proceed?

- 1. Take leaves of plants like peepal, rose, banana, bamboo, kaner, wheat, maize, mango, lemon and balsam.
- 2. Put each and every leaf separately between the folds of newspaper, and place thick heavy books over them.
- 3. Turn the leaves upside down daily, and if possible, change the newspaper, so as to avoid fungal infection. Take out all the leaves after a week. Now hold leaves one by one between your thumb and fore finger and study their venation. You can take the help of a magnifying glass, if needed. Fill in the table given below by writing the type of venation found in different leaves.

Sl.No.	Leaf	Type of venation
1.	Peepal	
2.	Rose	
3.	Banana	
4.	Bamboo	
5.	Kaner	
6.	Wheat	
7.	Maize	
8.	Mango	
9.	Lemon	
10.	Balsam	

GETTING TO KNOW PLANTS

. 35








Fig.6.7.2 Two types of leaf venation

- 1. The leaves have lines all over its lamina. These lines are called veins. The pattern of these veins differs in the leaves of different plants.
- 2. Some leaves have veins running parallel to each other (as in wheat/grass). Such a venation is called parallel venation. While others have interconnected veins forming a net like structure (like leaves of peepal, rose etc.,). This is called reticulate venation.

36 _____UPPER PRIMARY SCIENCE KIT MANUAL



Extension


1. Place the leaves of peepal in water for about a week . Take them out and put these in between the folds of newspapers and press them by putting some weight over them.

After a few days take them out. You can make your innovative cards by using these decorative leaves.

Notes			


GETTING TO KNOW PLANTS _____





How will you proceed?

- 1. Take a medium size and well watered plant.
- 2. Remove all the leaves of one twig of this plant.
- 3. Now take two polythene bags. One bag cover this twig and the other cover any other twig in which the leaves are intact.
- 4. Tie the mouths of the two bags around the bases of the two twigs.
- 5. Leave the plant in Sunlight for an hour or two.
- 6. Drops of water will appear on the inner side of the bag that covered the leafy twig.
- 7. No water drops will be seen on the inner side of the bag over the twig without leaves.



_____UPPER PRIMARY SCIENCE KIT MANUAL



- 1. The droplets of water appeared inside the polythene bag that covered the leafy twig.
- 2. These droplets are the condensed form of water vapours released during transpiration by the plants.
- 3. Plants release water into the atmosphere in the form of water vapours through their leaves. This process is called transpiration.

Extension

Repeat the above activity on a cloudy day.

Notes			

GETTING TO KNOW PLANTS _____







ACTIVITY 6.7.4

plant?

What is required? Why leaves are called Potted plant, black paper, the kitchen of the beaker, tripod stand, kerosene burner, iodine solution, water, petridish and methylated spirit.

How will you proceed?

- Destarch the leaves of a plant by keeping the plant in a dark room for about a few days.
- 2. Cover one of its leaves with black paper on which a design of a star is cut.
- Now place this plant in the sun. 3.
- Dip the leaf which was covered by black paper in boiling water for a minute to burst open the cells, to remove the entire chlorophyll.
- 5. Now transfer the leaf in a test tube with methylated spirit. Put this test tube in a beaker filled with water.
- Boil the water in the beaker over a burner. The leaf will become brittle and hard.




Fig.6.7.4 Sunlight is necessary for photosynthesis

After iodine test

_ UPPER PRIMARY SCIENCE KIT MANUAL



- 7. Place the leaf again in hot water to soften it.
- 8. Spread the leaf in a petri dish and pour iodine solution on it.

Observation

Only those parts of the leaf which were not covered by the paper, show the presence of starch. It is because only the cut out design could get the sunlight and prepared the starch in the presence of light.

What have you learnt?

The portions of leaf not covered with black paper synthesised food (starch) in the presence of sunlight due to which the colour of leaf turned blue, black in the presence of iodine.

Extension

Repeat the experiment with a variegated leaf.

GETTING TO KNOW PLANTS .

41





Activity 6.7.5

What are the major types of roots?

What is required?

Grass and small

herbaceous (mustard /

Amaranthus) plants.

How will you proceed?

- 1. Dig out a small herbaceous plant from ground.
- 2. Wash the soil of the roots and observe the roots.
- 3. What do these roots look like? Are they similar or different?




Fig.6.7.5 (a) *Tap root* (b) *fibrous root*

What have you learnt?

- 1. There are two major types of root system tap root and fibrous root.
- 2. Some plants have roots which consist of a main root with side branches (Tap roots).
- 3. There are some other plants which do not have any main root and all the roots originate directly from the base of the stem (Fibrous roots).

42 _____ UPPER PRIMARY SCIENCE KIT MANUAL





- 4. Amaranths has tap roots.
- 5. Grass has fibrous roots.

Extension

1. You have learnt that the leaves of some plants have a parallel venation, while others have a reticulate venation. Venation of leaves is correlated with the type of roots these plants have. Some common plants are given in the following table. Fill in the table given below by writing the type of venation and roots each plant possesses:

Sl.No.	Plants	Type of venation	Type of root
1.	Mooli		
2.	Gram		
3.	Meethi		
4.	Maize		
5.	Pudina		
6.	Bathuia		
7.	Dhania		


2. Visit your near-by garden and observe the plants and trees. Differentiate between herbs/shrubs and trees. Observe the type of venation in all these plants leaves.

GETTING TO KNOW PLANTS -

_ 43







How will you proceed?

- 1. Select two herbaceous plants of the same kind from ground and dig them out taking care that the roots do not break.
- 2. Put one plant in soil in one of the pot labelled as A.
- 3. Cut off the roots from the other plant and fix it in soil in another pot B.
- 4. Water them regularly and observe the plants after a week.




Fig.6.7.6. (a) Growing of the plant in soil with its root (b) plants in soil without roots

44 _____ UPPER PRIMARY SCIENCE KIT MANUAL



- Plant in pot A will continue to stand erect, and it is difficult to pull it out as the roots help in providing anchorage.
- 2. Plant in pot B will not stand erect and will fall and it is easier to pull it out.
- 3. The plant in pot A shows normal growth and development as roots also help in the absorption of water and minerals.
- 4. The plant in pot B will die as there are no roots.

Extension

Observe the roots of carrots, banyan tree, sugar cane and money plant, and comment on their other functions.

GETTING TO KNOW PLANTS

. 45





6.8 Motion and Measurement of Distances

ACTIVITY 6.8.1

How can we recognise periodic motion and circular motion?

What is required?

Pendulum bob, thread

(1 metre), wrist watch

or stop watch, G-clamp

and iron rod.

How will you proceed?

- 1. Fix the G-clamp at the edge of a table. Insert the iron rod in the hole of G-clamp such that about 10-15 cm length of the rod should be projected outside and tighten the screw [Fig. 6.8.1 (a)].

 G-Clamp
- 2. Tie a thread to the hook of the bob.
- 3. To make a pendulum, suspend the bob with the help of the thread from the rigid projecting iron rod.
- 4. When the bob is stationary, displace it a little from its initial position O to the new position A and then release.
- 5. Observe its motion.
- 6. What type of motion is it? Let us find out.
- 7. Note the time for 10 complete swings i.e motion along AOB and back along BOA of the pendulum.




Fig.6.8.1 (a) The motion of a simple pendulum

- 8. Repeat the step No. 7 three times.
- 9. Is the time taken for 10 swings the same in each case?
- 10. If the time taken is the same in each case, it means that pendulum repeats its motion after a fixed interval of time. So we can say that the motion of the pendulum is periodic.
- 11. Now move the bob in such a way that it moves in a horizontal circular path [Fig. 6.8.1(b)].




Fig.6.8.1 (b) The circular motion of the bob

- 12. What type of motion does the bob show now?
- 13. Should we call this motion as circular motion? Discuss.

- 1. Motion is said to be periodic if an object repeats its motion after a fixed interval of time.
- If a body moves in a circle, then it is said to be in circular motion.


Extension

- 1. Observe different motions around you and classify them as periodic and non-periodic.
- 2. Verify whether circular motion of the bob of the pendulum is periodic or not.
- 3. Also observe some objects making circular motions which are not periodic.

MOTION AND MEASUREMENT OF DISTANCES -







How can we measure the object using a thread?

What is required? Cotton thread, bangle, scale or measuring tape, pencil and

How will you proceed?

- Keep the bangle on a sheet of paper and draw its outer boundary.
- 2. Remove the bangle.
- 3. Take a point, say A, on the boundary drawn and mark it.
- Put one end of the thread at 'A' and fix it firmly with your thumb of the left hand.
- 5. Now place a small portion of the thread along the boundary keeping it taut with fingers of other hand to reach at point B as shown in [Fig.6.8.2].
- 6. Now hold the thread at point 'B' with one hand.
- 7. Using the other hand, taut a little more portion of the thread along the boundary.
- Go on repeating this process till you reach back to the point 'A'.




Fig.6.8.2 Measurement of the perimeter of the bangle

UPPER PRIMARY SCIENCE KIT MANUAL



- 9. Now make a mark on the thread where it touches again at 'A'.
- 10. Measure the length of the thread from one end to the marked point.
- 11. Repeat the steps from 4 10 at least one more time.
- 12. How is this measured length of the thread related to the perimeter of the bangle? This length gives the perimeter of the bangle.

We can measure the length of the curved objects by using a thread.

Extension

- 1. How can we measure the diameter of a circular object without knowing its centre? (Hint: The diameter of a circle is the largest distance between any two points on the boundary of the circle.)
- Use this method to find out the perimeter of rectangular and triangular objects.
 (a) Length of the curved line.

MOTION AND MEASUREMENT OF DISTANCES _

49





ACTIVITY 6.8.3

How can we recognise the rectilinear motion of an object?

What is required?

A glass marble, poster colour/ink, a long sheet of white paper and a measuring tape.

How will you proceed?

- 1. Spread the white sheet of paper on a levelled surface.
- 2. Dip the glass marble in ink or poster colour (thick solution).
- 3. Roll the coloured marble on the white sheet of paper and observe the path followed by the marble from the impressions of the tracing on the papersheet.
- 4. Do we observe the path of the marble straight throughout?
- 5. Could you identify the motion of the marble in this case? Is it rectilinear?




Fig.6.8.3 Motion of object along straight path

0 _____UPPER PRIMARY SCIENCE KIT MANUAL



6. Measure the distance travelled by the marble while it is in rectilinear motion.

What have you learnt?

When an object moves along a straight path its motion is called rectilinear.

Extension

Perform the above activity using a pingpong ball and a steel ball-bearing and observe in which case the motion is rectilinear for a longer time.

Notes			

MOTION AND MEASUREMENT OF DISTANCES ____



6.9 Light, Shadows and Reflections

ACTIVITY 6.9.1

What is required?

Do all substances

allow light to pass
through them?

paper, wooden scale, glass
sheet, glass beaker,
cardboard, white sheet of
paper and wash bottle.

How will you proceed?

1. Look at an object say a pencil through the plastic scale, tracing paper, wooden scale, glass sheet, glass beaker, cardboard, white sheet of paper and washbottle separately [Fig 6.9.1].




Fig.6.9.1 Seeing object using trasparent, translucent and opaque object



- 2. Group the objects through which the pencil
 - (a) is seen clearly
 - (b) can be seen faintly (dimly).
 - (c) is not seen at all.
- 3. Do you know the materials with which these objects are made of?
- 4. On the basis of above observations classify the materials as transparent, translucent and opaque.

some substances allow light to pass through them completely, some partially and some do not allow the light to pass through them at all. These are known as transparent, translucent and opaque substances respectively.

Extension


List some more examples of transparent, translucent and opaque materials which you find in your everyday life.

LIGHT, SHADOWS AND REFLECTIONS _

53







ACTIVITY 6.9.2


How can you make an opaque paper translucent?

What is required? sheet of opaque white paper, oil and an object(burning candle).

How will you proceed?

- Take a sheet of white paper.
- Look at an object say a burning candle through it. Can you see the candle flame Fig.6.9.2 (a)?
- Now smear few drops of oil at the middle portion of the sheet of paper.
- Look at the burning candle through this portion of the paper again.
- 5. What do you observe? Is the candle flame visible now [Fig.6.9.2 (b)]?
- 6. What changes did you notice to the portion of the white paper when you smeared oil on it? The portion of the paper becomes translucent and therefore you are able to see the object (candle

flame faintly) through it.






Fig.6.9.2 (a) and (b) Seeing object through a translucent paper prepared by smearing oil on it

What have you learnt?

A translucent paper can be prepared by smearing oil on a thin sheet of opaque white paper.

UPPER PRIMARY SCIENCE KIT MANUAL



ACTIVITY 6.9.3

How are shadows formed?

What is required?

Generator, a bob, torch

bulb, transparent glass

plate, a ground glass plate

and a metallic scale.

How will you proceed?

- 1. Place a torch bulb near a wall (or a screen) and connect it with the generator.
- 2. Place an opaque object such as a bob or a metallic scale in between the bulb and the wall.
- 3. Rotate the handle of the generator to give electrical power to the bulb so that the bulb is lighted.
- 4. What do you observe on the wall? [Fig.6.9.3].
- 5. Is there a shadow on the wall?
- 6. Now move the opaque objects away from the bulb.



Fig.6.9.3 Opaque objects cast shadows on a screen when light falls on them

LIGHT, SHADOWS AND REFLECTIONS _____

_ 55



- 7. What happens to the nature of the shadow formed?
- 8. Now replace the bob by a transparent glass plate.
- 9. Do you still observe a clear shadow on the wall?
- 10. Replace the transparent glass plate by a ground glass plate (translucent).
- 11. Could you see a clear shadow of the ground glass plate on the wall?
- 12. In which of the above three cases is the shadow most clearly formed?
- 13. On the basis of your observations, conclude the kinds of obstacles: (transparent/translucent/opaque) which can form a clear shadow.


which can form a clear shadow.	
What have you learnt?	
Opaque objects cast shadows on a screen more cl	'early
when light falls on them.	
Notes	

56 _____UPPER PRIMARY SCIENCE KIT MANUAL









How will you proceed?

- 1. Light the candle and keep it at about 2 m distance from you.
- 2. Keeping the tube straight in the direction of the candle flame, observe the flame from the other end of the tube [Fig 6.9.4 (a)].
- 3. Is the flame visible?
- 4. Then bend the tube in different directions and observe the flame again through the tube [Fig 6.9.4 (b)].
- 5. Is the flame visible now? Why is the flame not visible when the tube is bent?





Fig.6.9.4 Using straight and bent tubes to see objects

What have you learnt?

Light travels in a straight line.

LIGHT, SHADOWS AND REFLECTIONS ______5





ACTIVITY 6.9.5

How an image is formed in a pinhole camera?

What is required?

Two thick black paper sheets of size(15 cm X 20 cm), tracing paper/glue/cellotape, pair of scissors and a needle.

How will you proceed?

- 1. Make two cylindrical tubes using the thick black paper sheets such that the one just slides into the other.
- 2. Fix the tracing paper to one end of the inner sliding tube.
- 3. Cover one end of the outer tube with a circular piece of thick black paper sheet by using glue.
- Make a pinhole at the centre of the circular piece of black paper sheet attached at the end of the outer tube.
- 5. Slide the inner tube into the outer tube keeping the tracing paper towards you. Now, the pinhole camera is ready [Fig 6.9.5].
- 6. Keep the pinhole of the camera towards the well illuminated distant object like a tree or a building.
- 7. Observe how the distant object looks like on the tracing butter paper attached at the end of the inner tube.
- 8. Slide the inner tube inside the outer tube and find out the position when the image of the distant object is seen most clearly on the tracing paper.
- 9. Move the pinhole camera up, down and sideways while observing the image of the distant object.

UPPER PRIMARY SCIENCE KIT MANUAL





- 10. Do you still see the image of the distant object on the tracing paper of the pinhole camera in all the positions distinctly? If not, discuss why?
- 11. What do you conclude from the above observations?
- 12. Draw diagram to show how the pinhole camera works and state how light travels.

Inverted Image
Pinhole
Camera

Fig.6.9.5 Pinhole
camera


What have you learnt?

Light travels in a straight line.

LIGHT, SHADOWS AND REFLECTIONS __

59





- 1. Connect the bulb of the ray streak apparatus with the generator.
- 2. Place the ray streak apparatus on a white sheet of paper.
- 3. Rotate the handle of the generator.
- 4. What do you observe on the paper when the bulb starts glowing?
- 5. Is the beam of light parallel?
- 6. Place a plane mirror in its path [Fig 6.9.6].
- 7. What happens to the beam of light falling on the mirror?





Fig.6.9.6 A plane mirror reflects a beam of light

00 ______ UPPER PRIMARY SCIENCE KIT MANUAL



- 1. Light falling on a plane mirror is reflected.
- 2. A parallel beam of light falling on the plane mirror is reflected as a parallel beam.

Extension

Turn the plane mirror first towards left and then towards right and observe the reflected beam.

- 1. What do you conclude from the above observations?
- 2. Do you find any change in the shape of the reflected beam?
- 3. Do you find any change in the direction of the reflected beam?

LIGHT, SHADOWS AND REFLECTIONS -

61







How will you proceed?

- 1. Take a plane mirror.
- 2. Stand in front of it.
- 3. What do you observe in the plane mirror?
- 4. Raise your left arm and observe your image.
- 5. Which arm of the image appears to be raised left or right? Repeat this with your right arm and conclude [Fig. 6.9.7].



Fig.6.9.7 Image formation using plane mirror

What have you learnt?

- A plane mirror forms an image of an object placed in front of it.
- 2. In a mirror the left appears right and the right appears left such that only the sides are interchanged.

62 ______UPPER PRIMARY SCIENCE KIT MANUAL



ACTIVITY 6.9.8

What is the nature of the image formed by a plane mirror?

What is required?

Plane mirror, candle,

match box and a

mirror stand.

How will you proceed?

1. Place a plane mirror vertically on a table using a stand.

2. Now bring a lighted candle in front of the mirror. [Fig. 6.9.8].

3. Observe the image of the flame of the candle in the mirror.

4. Is the image formed in the mirror

- (a) erect or inverted?
- (b) enlarged or diminished?
- (c) of the same size as the object?
- 5. Now move the candle forward and backward and observe whether the image of the candle flame
 - (a) changes in size.
 - (b) remains erect or becomes inverted as the flame moves towards or away from the mirror.





Fig.6.9.8 Image formed by using a plane mirror

What have you learnt?

Image formed in a plane mirror is erect and is of the same size as the object.

LIGHT, SHADOWS AND REFLECTIONS _______







How can we see objects which cannot be seen directly?

What is required?
Periscope (from kit).

How will you proceed?

- 1. Take periscope box from the kit.
- 2. Keep one end of the periscope towards an object and see through the other end [Fig 6.9.9] by adjusting the angle of the mirror.
- 3. Is the object visible? If yes, why?
- 4. Now use the periscope to see objects placed on the table while you are sitting underneath and can not see them directly.
- 5. Explain why you can see the objects placed on the table while you are sitting underneath.



What have you learnt?

Fig.6.9.9 Periscope


Light can undergo multiple reflections.

Extension


Make your own periscope by using cardboard or two large tooth paste packings by fixing mirror strips.

64 ______UPPER PRIMARY SCIENCE KIT MANUAL





6.10 ELECTRICITY AND CIRCUITS



How will you proceed?

- Take a circuit board and connect the cell holder and torch holder in its.
- 2. Place the torch bulb and cell in the respective holder.
- 3. Also place the cell and switch in the sockets of the circuit board.
- Now connect the bulb to the cell through the switch as shown in [Fig.6.10.1]. This is a simple circuit having a dry cell, switch and a torch bulb.




Fig.6.10.1 A simple electric circuit

- 5. What happens to the bulb when circuit is switched ON?
- 6. What do you observe when circuit is switched OFF?

A simple electric circuit consists of a cell, bulb, a switch and connecting wires. When the switch is on, the circuit is closed and the bulb glows.

Extension

- Interchange the position of the bulb and the switch on the circuit board. Observe whether the bulb glows or not.
- Interchange the +ve and -ve terminals of the cell on the circuit board and observe its effect on the intensity of the glow of the bulb.
- 3. Note down your observations and conclude.

06 _____UPPER PRIMARY SCIENCE KIT MANUAL





ACTIVITY 6.10.2

How can we identify conductors and insulators?

What is required?

Kits items made up of different substances such as plastic scale, glass rod, eraser, iron rod, aluminium rod, piece of wood, rubber band, circuit board, banana pin plug, torch bulb and a dry cell etc.

How will you proceed?

- 1. Make a circuit on the circuit board with bulb, dry cell and switch as shown in [Fig.6.10.2].
- 2. Insert the objects (i.e plastic scale, glass rod, iron rod etc.,) one by one in the gap A B. Observe whether the bulb glows or not.




Fig. 6.10.2 Circuit board connected with a bulb, switch and a cell

ELECTRICITY AND CIRCUITS _____



3. Note down your observations and conclusions in the table given below.

Sl.No.	Name of the object	Material	Observation (whether the bulb glows or not)	Conclusion (Insulator / Conductor)
1.	Eraser	Rubber		
2.	Iron rod	Iron		
3.				
4.				

What have you learnt?


Bulb glows when current passes through it. Some materials allow electric current to pass through them and some do not. Conductors allow the electric current to pass through them and that is the reason why the bulb glows. Reverse is the case with insulators.

Extension

Repeat the above activity with a graphite rod and pencil lead. Does the bulb glow in each case? What do you infer from your observation?

68 ______ UPPER PRIMARY SCIENCE KIT MANUAL





ACTIVITY 6.10.3

How can we make an electric switch?

What is required? Wooden block or thermocol sheet, drawing pins and safety pin or paper clip.

How will you proceed?

- Insert a drawing pin into the ring of the safety pin.
- Fix it on the thermocol sheet or on the wooden block. 2.
- Make sure that the safety pin can be rotated freely.
- Now fix the other drawing pin on the thermocol sheet in such a way that the free end of the safety pin can touch it. The safety pin fixed in this way that it would serve as an electric switch.

5. Now take the circuit board having a cell and a torch bulb fixed on it.

With the help of connecting wire connect your electric switch in the circuit. [Fig. 6.10.3]

7. Observe the bulb when safety pin (a) touches the drawing pin (b) does not touch the drawing pin?

8. When does the bulb glow?



Fig. 6.10.3 Making an electric switch

What have you learnt?


An electric switch is a make and break arrangement. When contact is made we say switch is ON, the

circuit gets closed and the bulb glows.

ELECTRICITY AND CIRCUITS _

sheet





ACTIVITY 6.10.4

Does tap water conduct electricity?

What is required?
Two dry cells, cell
holder, one light
emitting diode(LED),
plug key, a plastic cup,
tap water, distilled
water connecting wires
and crocodile clip.

How will you proceed?

- 1. Fill tap water in a plastic cup.
- 2. Dip two crocodile clips connected with wires in water without touching each other.
- 3. Connect one of the wires joined with these clips with a LED and the other wire to the cell through a switch and complete the circuit as shown in the [Fig 6.10.4].
- 4. Now switch ON the circuit and check whether LED glows. If the LED does not glow, then interchange




Fig 6.10.4 Circuit arrangement using LED to test whether tap water conduct electricity or not

70 _____ UPPER PRIMARY SCIENCE KIT MANUAL



- the terminals of the LED and check again [Note: LED glows only when current flows through it in a particular direction]. Does the LED glow now? What does this indicate?
- Does the tap water conduct electricity or not? Give your conclusion.

Tap water conducts electricity.

Extension

Repeat the above activity with distilled water and observe. What do you infer? Add some salt to this water and repeat the activity. What inference can you draw now?

Notes			

ELECTRICITY AND CIRCUITS -



6.11 Fun with Magnets

ACTIVITY 6.11.1

How can we test
whether a given object
is a magnet or not?


What is required?

Magnets of different
shapes, a pencil, an
iron nail, a glass rod,
thread, iron filings, Gclamp alpins, a sheet
of paper, shoe nails
and a compass needle.

How will you proceed?


1. Place the objects (i.e. pencil, iron nail, glass rod, magnets, shoe nail, alpins etc.), one by one on iron filings spread on a paper.

2. See if the iron filings stick on the different objects or





- 3. If the iron filings stick to the object, it is a magnet, otherwise not.
- 4. Classify the above objects in two groups:(a) Magnets and (b) Non-magnets
- 5. Alternatively, suspend these objects one by one from a rigid support using a thread.
- 6. Observe in which direction the object rests.
- 7. Now disturb it from its position of rest and wait till it again comes to rest.
- 8. Note the direction in which it settles again.
- 9. Does the object settle in the same direction as previously? If yes, it is a magnet. If not, it is not a non-magnet.



- 1. Magnet attracts iron.
- 2. A freely suspended magnet always rests in a fixed direction.

FUN WITH MAGNETS _____

_ 73







Extension


Bring each object one by one near to one end of a compass needle. Note whether the object attracts or repels that end of the needle.

Now bring the other end of the object near to the same end of the compass needle and note the difference. What conclusion can you draw from this activity in the light of the fact that compass needle is also a small magnet and similar poles of two magnets repel each other?

Notes				

74 ______UPPER PRIMARY SCIENCE KIT MANUAL





What is required?

Magnets of different shapes, a sheet of white paper and a pencil.

How will you proceed?

- 1. Take all the magnets given in the kit.
- 2. Arrange all the magnets on a table.
- 3. Observe the shapes of these magnets. Do they look alike? If no, then arrange them in different groups according to their shape [Fig. 6.11.2].
- 4. From your observation, how many different shapes of magnets do you have now on the table?
- 5. Observe the shape of magnets in each group.



Fig 6.11.2 Different types of magnets

FUN WITH MAGNETS _____

_ 75



- 6. Draw a rough diagram of each of these shapes on a sheet of paper.
- 7. Suggest a name for the shape of magnets in each group.
- 8. Record your observations in the following table.

Group	Shape of the magnets as they appear	Suggested name of the magnet according to its shape
1.		
2.		
3.		
4.		
5.		

Magnets can have different shapes.


Extension

Observe different devices/ toys in your every day life. Try to find out if there is any magnet used in them. Make a list of devices/toys where magnets are used. Also try to draw the shapes of these magnets.

76 ______UPPER PRIMARY SCIENCE KIT MANUAL







ACTIVITY 6.11.3

Do all materials experience mutual attraction with a magnet?

What is required?

A strong magnet, objects made of different materials such as wood, plastic, iron, aluminum, glass, rubber and nickel etc.

How will you proceed?

- 1. Place the given objects made of different materials on a table.
- 2. Now bring a strong magnet one by one near each of the object. What do you observe?
- 3. Do you feel any attraction between them in each case?
- 4. Do you find that the magnet and the objects tend to move towards or away from each other?
- 5. Do you find that there is no force between the magnet and some of the objects at all?




Fig.6.11.3 Mutual attraction of a magnet using different materials

FUN WITH MAGNETS _______77



Record your observations in a tabular form. A sample table is given below for reference.

Material of	Observation	Conclusion
the object	(attraction/	(magnetic substance/
	repulsion/	non-magnetic
	neutral)	substance)
Plastic		
Rubber		
Iron		

What have you learnt?						
Materials attracted by	a magnet are known as					
magnetic materials.						

Notes			

78 _____ UPPER PRIMARY SCIENCE KIT MANUAL





Does a magnet attract iron filings / nails equally at all its points?

What is required? A bar magnet, U-shaped magnet, iron filings and a nails or alpins etc.

Wrap the magnet in paper so that iron filings can be easily removed afterwards.

How will you proceed?

- Spread iron filings on a sheet of paper placed over a table.
- 2. Put a bar magnet over the filings.
- Take out the magnet. Observe the iron filings clinging to the magnet on its surface.
- Do you find the iron filings uniformly stuck on the surfaces of the magnet? If not, where do you find the iron filings sticking in larger quantity?
- Draw a diagram to represent your observation.
- Compare your diagram with [Fig 6.11.4 (a)].




Fig.6.11.4 (a) Iron filing sticking to a bar magnet

FUN WITH MAGNETS _



- 7. Repeat this activity by using a U-shaped magnet. What is your observation now?
- 8. Again draw a diagram to show your observations using U-shaped magnet.
- 9. Now, bring few iron nails near different points of one end of the bar magnet so that they stick one behind the other as shown in [Fig. 6.11.4 (b)]
- 10. How many nails can you stick at different points on the bar magnet like that ?
- 11. Can you stick the same number of nails at every point of the magnet? Observe and give your conclusion.




Fig.6.11.4(b) More nails stick to the magnet near its end

The amount of iron filings clinging is maximum at the two ends of a magnet.

0 _____UPPER PRIMARY SCIENCE KIT MANUAL



ACTIVITY 6.11.5


How do we names the poles of a magnet?

What is required?

Compass needle, bar magnet, U-shaped magnet, a thread, a piece of chalk and G-clamp.

How will you proceed?

- 1. Draw a straight line on a table or floor using a chalk along north-south direction.
- 2. Place a compass needle over this line.
- 3. Observe the direction along with the needle rests.
- 4. Slowly rotate the compass keeping the centre of its needle on the line.
- 5. Observe the direction of the needle in different positions.
- 6. Does the needle of the compass change its alignment?
- 7. Can you make it align in a different direction?
- 8. Which end of the compass needle is pointing towards north?






Fig.6.11.5 Finding the poles of a magnet using compass needle

FUN WITH MAGNETS _

_ 8′



- 9 Which end of the compass needle is pointing towards south?
- 10. Does the rotation of the compass change the direction of the needle as shown by these ends?
- 11. Now suspend a bar magnet freely from a stand using a thread.
- 12. Observe the direction in which the magnet rests.
- 13. Which end of the magnet points towards north?
- 14. Which end points towards south?
- 15. Mark 'N' at the end of the magnet seeking towards north. This is the 'North Pole' of the magnet.
- 16. Write 'S' at the end of the magnet seeking towards south. This is the 'South Pole' of the magnet.

The end of a freely suspended magnet pointing towards north is north pole and pointing towards south is south pole.


Extension

Repeat the above activity using a U-shaped magnet.

_____UPPER PRIMARY SCIENCE KIT MANUAL







ACTIVITY 6.11.6

Can we make a magnet by ourselves?

What is required?

Bar magnet, iron

nail, sewing needle

and iron strip.

How will you proceed?

- 1. Take an iron nail and dip it in iron filings.
- 2. Check whether it is already a magnet or not.
- 3. Make sure that the nail chosen is not already a magnet .
- 4. Now place the nail on a flat surface.
- 5. Rub the nail with a bar magnet by touching its one pole along its length [Fig.6.11.6].




Fig 6.11.6 Rubbing a mgnet over iron nail

- 6. Rub the nail for some time.
- 7. Now test whether the nail has become a magnet or not.
- 8. You can test this by dipping the nail in iron filings.
- 9. What is your observation?
- 10. Does the nail become a magnet now?

What have you learnt?

An iron nail can be magnetised by rubbing a magnet on its surface.

Extension

Repeat the above activity to magnetise other iron objects like an iron strip or a sewing needle.

FUN WITH MAGNETS ____







How can we make a simple compass needle?

What is required?

Bar magnet, sewing needle, thermocol disc of about 4 cm diameter and 1 cm thickness, small plastic tub, water and adhesive tape.

How will you proceed?

- 1. Magnetise the needle by rubbing it with a pole of a strong bar magnet (as suggested in the previous activity).
- 2. Place the needle at the centre of the thermocol disc and fix it with adhesive tape.
- 3. Place it gently over water taken in a tub such that it floats [Fig 6.11.7].
- 4. Observe the direction along which the needle rests. Is it pointing towards north-south direction?
- 5. Now your compass is ready. Sewing needle





Fig.6.11.7 Making a simple compass needle

What have you learnt?

A needle can be magnetised and can be used as a compass needle.

84 ______UPPER PRIMARY SCIENCE KIT MANUAL





- 1. Take two bar magnets.
- 2. Hold one magnet in each hand.
- 3. Bring the north-pole of one of these magnets near the north pole of the other magnet.
- 4. Observe what happens.
- 5. Now bring the south-pole of one of these magnets near the south-pole of the other.
- 6. What happens now?
- 7. What do the magnetic poles experience?
- 8. How do they behave? Does the similar pole of two magnets attract or repel each other?
- 9. Now, bring the north-pole of one of these magnets near the south pole of the other.
- 10. Note the change in the type of force experienced by the two magnet. Do you observe the force of attraction or repulsion between them?
- 11. Tabulate your observations.




Fig. 6.11.8 (a) Unlike poles attract each other (b) Like poles repel each other

FUN WITH MAGNETS _____

What	have	vou l	learnt?
rriuu	mu e	you i	em ru.

Like poles repel and unlike poles attract each other.

Extension

Repeat the above activity using the cylindrical magnets. What do you observe? Is it different from your observations with respect to the bar magnets?

Notes			

86 _____UPPER PRIMARY SCIENCE KIT MANUAL



ACTIVITY 6.11.9

How can we compare the strength of two magnets using a compass needle?

What is required?
Two bar magnets and
a compass needle.

How will you proceed?

- 1. Place a compass needle on a table.
- 2. Mark its boundary.
- 3. Notice that the compass needle points in north-south direction.
- 4. Now place the north pole of a bar magnet at a certain distance from the compass needle.
- 5. Note the extent of deflection on the compass needle.
- 6. Remove the bar magnet away from the compass needle.
- 7. Bring another bar magnet towards the compass needle.
- 8. Place its north pole exactly in the same position and at the same distance as taken in step
 No. 4 from the compass needle.
- 9. Again note the deflection on the compass needle.
- 10. Which magnet causes greater deflection on the compass needle?

Compass needle (a) Lower deflection



Magnet causing more deflection on the compass

needle is having

more strength.



Fig. 6.11.8 Finding the strength of magnet using compass needle

FUN WITH MAGNETS _____





ACTIVITY 6.11.10

Can a magnet lose its
magnetic properties
when heated?

What is required?

A bar magnet,
forceps, kerosene
burner, iron filings,
compass needle and a
sewing needle.

How will you proceed?

- 1. Take a sewing needle.
- 2. Rub it with a pole of a bar magnet for some time.
- 3. The needle becomes a magnet now.
- 4. Verify this with the help of a compass needle.
- 5. Also, verify this by dipping the needle in iron filings.
- 6. Make a rough estimate of the amount of iron filings that cling to the needle.
- 7. Remove all the iron filings from the needle.
- 8. Hold the needle over a flame with the help of forceps. Heat it for about 8 -10 minutes.




Fig.6.11.10 (i) Rubbing magnet over sewing needle

88 _____UPPER PRIMARY SCIENCE KIT MANUAL






Fig.6.11.10 (ii) Magnet loses its strength on heating

- 9. Allow the needle to cool.
- 10. Now dip it again in iron filings.
- 11. Make an estimate of the amount of iron filings that cling to it now.
- 12. What is your observation?
- 13. What conclusion can you draw from this?
- 14. Can you test the strength of the needle-magnet using a compass needle?


Magnet loses its magnetic property on heating.

Extension


Check whether a needle magnet also loses its magnetic property by hammering or not.

FUN WITH MAGNETS _____





6.12 WATER



- 1. Take two clean watch glasses and label them as A and B.
- 2. Keep A under direct sunlight and B in the open but under the shade.
- 3. Put 1 mL of water in both the watch glasses with the help of a syringe.
- 4. Occasionally observe both of them.





Fig.6.12.1 Evaporation in sunlight is faster than in shadow

90 _____ UPPER PRIMARY SCIENCE KIT MANUAL



- 5. Did you find any difference in the amount of water kept in them?
- 6. From which plate, did the water disappear faster?
- 7. Now touch both the watch glasses.
- 8. Which one was warmer? Can you think of the reasons?
- 9. Also, did you notice any mark or residue left in the watch glass after complete evaporation? What is this due to? Did water contain some dissolved materials?

- 1. Evaporation is the process of conversion of a liquid (water) into its vapour (water vapour). It depends on various environmental factors.
- 2. Water kept in the watch glass A, which was kept under direct sunlight disappears faster, as compared to that in watch glass B kept under shade. Watch glass A becomes warmer than watch glass B.
- 3. The process of evaporation becomes faster at higher temperatures. In the above activity, water kept in watch glass A received direct and strong heat from the sun. This heat helped in water getting quickly evaporated. On the other hand, water in watch glass B was kept in the shade could get heat only from the warm air. Thus, though water evaporated from this watch glass also, but the process was much slower.

WATER _____







Extension

- 1. Repeat the above activity by keeping watch glass A in a closed room under circulating fan and watch glass B in another closed room without fan. Observe carefully from which watch glass does water disappear faster? Give reasons for your observations.
- 2. Take petri dishes which differ in size. One may be of 5 cm diameter and another of 7 cm diameter. Take equal volumes of water (say 5 mL) in both. Keep both of them under direct sun light. Observe the disappearance of water. In which dish does water disappear faster?

Related Facts


- 1. Wet clothes dry sooner under the direct heat of sun as compared to those kept under shade.
- 2. In a closed room, wet clothes dry faster under a fan as compared to those kept in a room in which there is no circulation of air.
- 3. Water evaporated faster if taken in a wider vessel.
- 4. Rivers or lakes sometimes dry up during summers.
- Freshly mopped floors dry very quickly in summers than in winter and under a running fan than in a room with a switched OFF fan.

92 ______UPPER PRIMARY SCIENCE KIT MANUAL









Air contains water vapours which enter the air through the processes of **evaporation** and **transpiration**.

How will you proceed?

- 1. Take a beaker (clean its outer surface with a piece of cloth if necessary) and keep it on a flat surface.
- 2. Put three or four ice cubes in the beaker.
- 3. Observe the outer surface of the beaker. What do you notice?
- 4. You must have noticed that the water droplets appear on the outer surface of the glass.
- 5. Where has this water come from?




Fig.6.12.2 Process of evaporation and transpiration

WATER ______93



- 1. The air contains water vapours.
- 2. When the water vapours come in contact with cooler outer surface of the glass, it condenses on it in the form of water droplets.

Related facts

- 1. As we go higher up in the atmosphere, the air becomes cooler. Because of this, the water vapours present in the air condense to form clouds.
- 2. When we take out any object from the refrigerator, water droplets are formed on its outer surface after some time.

Notes			

94 _____UPPER PRIMARY SCIENCE KIT MANUAL





6.13 AIR AROUND US

ACTIVITY 6.13.1

What can the moving air do?

What is required?

A square piece of paper of size (22 cm × 22 cm), a pencil, a pair of scissors, a drawing pin, a small stick, crayon colours.

How will you proceed?

A. How to make the firki?

- 1. Take a square piece of paper of size ($22 \text{ cm} \times 22 \text{ cm}$).
- 2. Fold the paper in half twice as shown in the Fig. (a) and (b) and you get a small square.
- 3. Fold the square into triangles as shown in the Fig. (c).
- 4. With the help of scissors, cut the middle line of the triangle upto three fourth so that the middle portion does not get separated as shown in Fig. (d).
- 5. Open the folds and put the numbers 1 8 on the triangles or put some colours on them as shown in Fig. (e).
- 6. Then fold the alternate number tips to the centre as shown in Fig. (f).
- 7. Fold 1, 3, 5, 7 to the centre and hold them together as shown in Fig. (g).
- 8. Fix these folded tips of paper to a stick by using a drawing pin and cover the pointed end of the pin by using plastic so that it does not hurt any one.





Go to an open place and hold the stick. Turn the firki in different directions and move it a little back and forth and observe what happens?




Fig. 6.13.1 Arrangement to make the firki using square sheet of paper



What have you learnt?

The moving air helps in rotating the firki.

Extension


You can try this activity by taking the firki in front of a table fan and also in a room without any or very small opening for the entry of air.

Related facts

- 1. Moving air rotates the blades of a wind mill.
- Air helps in the movement of sailing yachts, gliders, parachutes etc.
- Very strong winds can cause cyclones, storms etc. and which can have destructive effects.

_____UPPER PRIMARY SCIENCE KIT MANUAL





Air is present everywhere around us. It occupies all the space. It does not have any colour and is transparent.

How will you proceed?

- 1. Take a bucket and fill it about 3/4th with water.
- 2. Now take an empty plastic bottle and hold it upside down.
- 3. Try to fill the bottle by dipping it into the water taking care that it should be absolutely vertical.
- 4. What did you observe? Did water enter the bottle?
- 5. Now, tilt the bottle a little and try to fill the water in the bottle. What do you observe?
- 6. Record your observations and think of the reasons.




Fig.6.13.2 Presence of Air in empty bottle

What have you learnt?

Air is present everywhere and occupies all the available space. The bottle which looks empty is actually full of air. When the bottle was dipped in

AIR AROUND US _____



an inverted and vertical position, water did not enter the bottle because air could not escape from it. But when the bottle was slightly tilted, air escaped from the bottle which was observed in the form of bubbles. Also, it created some space in the bottle which was occupied by the water and the bottle got filled with water.

Extension

- 1. Place a ping-pong ball on the surface of water in a bucket. Push the ball downwards with the help of the mouth of an inverted bottle. Observe carefully. Now slightly tilt the bottle and observe.
- 2. Take a bottle and fill it about half with water. Now, take a small plastic sheet and cut in the form of diver. Place this diver on the surface of water. Cover the mouth of bottle with a rubber sheet. Press the rubber sheet and observe. Give reasons for your observations.
- 3. Find out how a hovercraft moves.

Related facts

- 1. All the vehicles move on tyres which are filled with the compressed air.
- 2. Mountaineers or divers carry cylinders with them which are filled with oxygen.

98 _____UPPER PRIMARY SCIENCE KIT MANUAL





ACTIVITY 6.13.3

What does air contain?

Air is a mixture of gases. Of all these gases, burning can occur only in the presence of oxygen. What is required?

A plastic or glass tub, two small candles of same size, two glass tumblers of different sizes but both larger than candles, some water, a match box.

How will you proceed?

- 1. Take the glass or plastic tub. Fix two candles in it leaving a distance of approximately 6 cm in between them.
- 2. Now, fill the tub with some water.
- 3. Light the candles with the help of a matchstick.
- 4. Now take two glasses of different sizes and invert them on each of the candles. (Take care that candle height should be such that it should not go beyond the middle of the glass) [Fig.6.13.3].
- 5. Observe the burning candles and the water levels in each glass carefully.
- 6. Which of the candles stopped burning earlier and why?
- 7. Did you notice any rise in water level in the glasses? Was it same or different in the two glasses?
- You would have noticed that the flame of the candle goes out in the smaller glass tumbler earlier than that of the candle kept under bigger glass. Also, the

AIR AROUND US ______9





level of water rises more in bigger glass tumbler than the smaller glass tumbler. Can you think of the reasons?




Fig. 6.13.3 Arrangement to show air is a mixture of gases

What have you learnt?

Air is a mixture of gases. It mainly consists of oxygen(21%), nitrogen (78 percent), water vapour, carbon dioxide and other gases. When the burning candle is covered by an inverted glass tumbler, it uses up one component of air, i.e. oxygen. When most of the oxygen is used up by the candle, it stops burning and the flame goes out. The consumption of oxygen creates space within the glass tumbler, which gets filled by the water present in the trough. Water cannot fill the tumbler completely because rest of the gases present in air occupy the space inside the glasstumbler.

100 _____ UPPER PRIMARY SCIENCE KIT MANUAL



Extension

The above activity can be performed by putting a little red ink into the water.

This will make the visualisation of rising water level more clear.

Note for Teachers

The teacher may avoid discussing the above activity based on the per cent of oxygen present in the air.

The candles may not utilise 100 per cent oxygen and may blow out. As a result, the rise in water level may vary depending upon the conditions.

Related Facts

We inhale atmospheric air which contains oxygen. This oxygen reaches the cells and combines with the food (glucose). The oxygen 'burns' the food and releases energy.

AIR AROUND US _______10







UPPER PRIMARY SCIENCE KIT MANUAL

FOR CLASS VII






राष्ट्रीय शैक्षिक अनुसंधान और प्रशिक्षण परिषद् NATIONAL COUNCIL OF EDUCATIONAL RESEARCH AND TRAINING





7.1 Nutrition in Plants




tubes, tripod, water, rectified spirit, kerosene burner, petri dish, iodine solution, dropper, black paper, U-clip to hold the paper and potassium hydroxide solution

How will you proceed?

To test a leaf for starch (Iodine test).

- 1. Dip the leaf in boiling water for a minute to rupture the cells.
- 2. Boil the leaf in methylated spirit over a water bath (beaker filled with water) till it becomes pale-white due to the removal of chlorophyll. The leaf now becomes hard and brittle.
- 3. Place it again in hot water to soften it.






Before iodine test



After iodine test

Fig.7.1.1(b) Cholrophyll is necessary for photosynthesis



Before iodine test



After iodine test

Fig.7.1.1 (c) Only the parts of the leaf exposed to sunlight give the positive test with iodine test

- 4. Spread the leaf in a petridish and put 5 8 drops of iodine solution on it. Presence of starch will be indicated by a blue-black colour. A part of the leaf without starch will show brown coloration.
- **A** To show that chlorophyll is necessary for photosynthesis:
 - 1. Take a plant with (coloured) leaves having some green and some coloured areas. Example: Coleus and Croton.
 - 2. Destarch the leaves by keeping the plant in a dark room for a few days.
 - 3. Place the plant in the sun for a few hours.
 - 4. Afterwards, pluck one leaf.
 - 5. Make its outline on paper and mark inside the outline the green and coloured areas.
 - 6. Test the leaf for starch .Only the green parts of the leaf turn blue-black, showing the presence of starch.
- **B** To show that sunlight is necessary for photosynthesis.
 - 1. Take a plant with destarched leaves.
 - 2. Cover one of its leaves with black paper on which design is cut.
 - 3. Take the plant in the Sun for few hours.

UPPER PRIMARY SCIENCE KIT MANUAL



4. After a few hours, test the leaf which is covered by black paper for the presence of starch. It will be observed that only those parts of the leaf which were left uncovered by the paper, turn blue-black, showing the presence of starch.




Fig.7.1.1 (d) Carbon dioxide is necessary for photosynthesis

- **C** To show that carbon dioxide is necessary for photosynthesis:
 - 1. Take a branch of a plant with destarched leaves.
 - 2. Insert one of its leaves (through a split cork) into a conical flask containing potassium hydroxide. Potassium hydroxide absorbs carbon dioxide.
 - 3. Leave the plant in sunlight.
 - 4. After a few hours, test this and any other leaf of this plant for starch by using the iodine test.

 The leaf which was exposed to the atmospheric air becomes blue-black, and the one inside the flask containing potassium hydroxide does not become blue-black, showing that carbon dioxide is necessary for photosynthesis.

What have you learnt?

Plants manufacture their food (starch) with the help of water and carbon dioxide. This process (photosynthesis) occurs in the chlorophyll containing cells of the leaves and takes place in the presence of sunlight.


NUTIRITION IN PLANTS _

105





7.2 HEAT



How will you proceed?

1. Take water in a beaker.

- 2. Fix the thermometer in a G-clamp stand.
- 3. Look at the markings of thermometer, i.e, note the minimum and maximum temperature of the thermometer scale.
- 4. Dip the bulb of the thermometer in water kept in a beaker. Ensure that the bulb does not touch the wall of the beaker.




Fig.7.2.1 Arrangement for measuring temperature of water using a laboratory thermometer



- 5. Note the reading of the thermometer when mercury thread becomes steady. Remember that you should read the thermometer keeping the level of mercury along the line of sight.
- 6. Now heat the water with the help of the kerosene burner.
- 7. Put the beaker with hot water on the table. Now dip the thermometer in water and observe the mercury thread [Fig. 7.2.1].
- 8. Observe the mercury thread for about 10 minutes. What do you find?
- 9. Do you find any change in the reading of the thermometer?

Temperature of an object can be measured with the help of a laboratory thermometer.

Extension


Try to measure the body temperature using a laboratory thermometer.


Do you find any difficulty?

HEAT ________107









How will you proceed?

- 1. Wash the thermometer with clean water and then with antiseptic solution.
- 2. Note the reading of thermometer. It should be below 35° C.
- 3. If it is not below 35°C, then hold the thermometer firmly and shake it till the thread of the thermometer falls below 35°C.
- 4. Place the bulb of the thermometer under your tongue and keep it there for about one minute.
- 5. Now take out the thermometer from your mouth.
- 6. Note the reading of the thermometer in degrees Celcius by keeping the level of mercury along the line of sight.




Fig.7.2.2 *Measuring body temperature using clinical thermometer*

08 _____ UPPER PRIMARY SCIENCE KIT MANUAL



We can measure the body temperature by clinical thermometer.

Extension

- 1. Measure the body temperature of some of your friends in °C and then find the average body temperature. Compare this temperature with your body temperature.
- 2. Also, measure your body temperature in ${}^{\circ}F$.

Notes			

HEAT ________109



ACTIVITY 7.2.3

How can we show that the rate of heat conducted by different substances is different?

What is required?

T-shape holder, two rods of same dimensions made of different materials (such as iron and aluminium), kerosene burner, match box and candle wax, clamp, iron rod, boss head, G-clamp.

How will you proceed?

- 1. Take two identical rods of same dimension, one of aluminium and the other of iron.
- 2. Fix them in T-shape holder. Fix the holder in the clamp.
- 3. Fix small pieces of candle wax on each rod at equal distances as shown in the Fig.7.2.3.
- 4. Now heat the middle part of the T-shape holder.
- 5. What happens to the small pieces of wax fixed on the rods?




Fig.7.2.3 Arrangement to show the conduction of heat in aluminium and iron rod fitted with a T-shape holder



- 6. Observe from which rod does the wax start falling first?
- 7. How would you conclude your observations?

- 1. The rate at which heat is conducted by different substances is different.
- 2. Aluminium rod conducts heat more easily than an iron rod of the same dimensions.

Extension

- Observe the cooking utensils. Which metals are these made up of? Make a list of them.
- 2. Observe the handles of different cooking utensils. Name the materials used in making them.

HEAT 11



ACTIVITY 7.2.4

Which of the two identical bodies (one dark-coloured and the other light-coloured) will absorb more amount of heat when kept under similar conditions?

How will you proceed?

- 1. Paste a black sheet of paper on one can and a white sheet of paper on another can. Ensure that there is no gap in between paper and the tin can.
- 2. Pour equal amount of water in both the cans.
- What is required? Two laboratory thermometers, black sheet of paper, white sheet of paper, two identical tin cans, cellotape and water.
- In the beginning the temperature of water in both the cans must be same. To verify it, you can use laboratory thermometer.
- 4. Cover the cans from the top and leave both the cans in Sun for about two hours. Ensure that the two cans




Fig.7.2.4 Measuring of temperature



- must be covered by two pieces of small cardboard so that it will not obstruct in heating the cans.
- 5. Measure the temperature of water in both the cans.
- 6. Do you observe any difference in temperature of water in two cans?
- 7. What do you infer from this?

Dark coloured bodies are good absorbers of heat.

Extension

- Take two identical cloths/papers, one black and the other white. Leave them in the Sun for about 2 hours. Measure the temperature of the cloths/papers with the help of laboratory thermometer.
 - Do you find any difference in the temperature between the two?
- 2. Take two identical thermometers. Wrap the bulb of one thermometer with black cloth and the bulb of the other thermometer with white cloth. Leave both the thermometers in the Sun for about 20 30 minutes.


Do you find any difference in the readings of two thermometers?


What conclusion can be drawn from the above activity?

HEAT _______113









How will you proceed?

- 1. Take a beaker.
- 2. Fill it 2/3rd with water and place it on the tripod stand.
- 3. Place a crystal of potassium permanganate gently at the bottom of the beaker using a straw.
- 4. Now heat the water by keeping the kerosene burner just below the potassium permanganate crystal. Ensure that the flame of the burner is kept low so that heating is slow.




Fig.7.2.5 Transfer of heat in liquid



- 5. Note the changes taking place in the beaker carefully.
- 6. Do you observe any convection currents in the beaker?

In liquids heat is transferred by convection.

Extension

Perform this activity by using used tea leaves in place of potassium permanganate crystals. You may observe the movement of tea leaves when water is heated.

Notes			

HEAT _______115



Activity 7.2.6

How can we show heat transfer in gases?

What is required?

A sheet of paper, a match box, a candle, thread, a pencil and a pair of scissors.

How will you proceed?

- 1. Cut a round sheet of paper of diameter about 20 cm.
- 2. Draw a spiral on it as shown in Fig.7.2.6.
- 3. Cut the paper along the marked line using a pair of scissors.
- Suspend the paper from its central point vertically using a thread.
- 5. Now light a candle and bring it below the paper spiral. Ensure that the flame should




Fig.7.2.6 Arrangement to show heat transfer in gases

be kept at a distance so that the paper will not burn.

- 6. Observe what happens?
- 7. Will the paper spiral start rotating?
- 8. Why does the paper spiral start rotating?


What have you learnt?

Transfer of heat in gases takes place by convection.

The air near the candle become, light on heating and starts rising. The rising air strikes the paper spiral which starts rotating.




7.3 Acids, Bases and Salts



Substances can be acidic, basic or neutral in nature. Their nature can be identified using indicators. The substances which change their colour when added to acidic or basic substances are called indicators. Some common indicators are litmus paper, phenolphthalein, turmeric and extract of china rose flower.

How will you proceed?

- Dilute one drop of shampoo with 10 drops of water.
- Take red and blue litmus papers and turmeric paper.
- 3. Using a dropper, put one drop of shampoo solution on each of these indicator papers [Fig.7.3.1 (a)].
- What do you observe? Record your observations in the table.










Fig.7.3.1 Test the acidic, basic or neutralnature of a substance

- 5. Put 5 drops of shampoo solution into a test tube containing some china rose extract [Fig.7.3.1 (b)].
- 6. What change in colour do you observe? [Fig.7.3.1 (b)]. Record your observations in the table given below.

Sl. No	Test Solution		Effec	t on	
		Red litmus paper	Blue litmus paper	Turmeric paper	China rose extract
1.	Shampoo				
	solution				
2.					
3.					



Shampoo is basic in nature because it turns

_ Red litmus blue.
_ Turmeric paper red.
_ China rose indicator green.

Extension


- 1. Test the nature of some more substances like soda water, lemon juice, soap solution, vinegar, sugar solution, common salt, dil. HCl, dil. H₂SO₄. It sodium hydroxide, ammonium hydroxide using the above indicators. Record your observations in the given table.
- 2. Prepare your own natural indicator using red cabbage leaves. Boil the cabbage leaves and strain them to obtain your own natural indicator. What colour does your indicator give with acids and bases?

ACIDS, BASES AND SALTS

119







How will you proceed?

- 1. Take 1 mL of dilute HCl in a test tube using a graduated dropper.
- 2. Add one drop of phenolphthalein into it and shake the solution. Do you observe any change in the colour of the solution?
- 3. To the above solution add 1-2 drops of sodium hydroxide solution by a dropper and shake the tube gently. Is there any change in the colour of the solution?
- 4. Continue adding sodium hydroxide solution drop wise with stirring till you observe the change in colour. Which colour is obtained?




Fig.7.3.2 Process of neutralisation



- 5. Now add 1-2 drops of dil. hydrochloric acid. Do you again observe any change in colour?
- 6. Again add 1-2 drops of sodium hydroxide solution. What happens to the colour of your solution?

When an acidic solution is mixed with a basic solution both the solutions neutralise each other as indicated by the change in the colour of the indicator. This reaction is known as neutralisation reaction.

Acid + Base \rightarrow Salt + Water

Note that in this reaction, salt and water are obtained as products

Extension

Observe the colour of the acid and base solutions before neutralisation and colour of the neutralised solution after mixing them using red and blue litmus papers.

ACIDS, BASES AND SALTS

121







7.4 Physical and Chemical Changes

ACTIVITY 7.4.1

How will you get water from ice and ice from water? How will you get water from water vapour?

How will you proceed?

- 1. Take some ice in a beaker and place it on the tripod stand and heat it with the help of kerosene burner.
- 2. What do you see? Ice melts to form water.

What is required?
Ice, beaker,
kerosene burner,
watch glass.




Fig.7.4.1 Ice melts into water on heating and water gets evaporated to form water vapour



- 3. Can we use this water to get back ice?
- 4. Suggest a method for it. Can you say that it is a physical change?
- 5. Take some water in a beaker. Heat the beaker over the kerosene burner. What do you see?
- 6. Keep a watch glass inverted above the beaker.
- 7. Do you see some water droplets on the inner surface of the watch glass?
- 8. You will see the droplets falling back into the beaker.
- 9. What is this phenomenon called as?

Conclusion

Ice melts to water and water can be frozen back to ice. Water gets evaporated to form water vapours which can condense to form water.

Extension


Think of some way to prove presence of water vapours in air.

PHYSICAL AND CHEMICAL CHANGES _

123







Sometimes a physical change involves change in colour.

How will you proceed?

- 1. Hold a shaving blade with a pair of tongs over the flame of a kerosene burner.
- 2. Does the colour of the blade change?
- 3. Remove the blade from the flame.
- 4. What do you observe after sometime?

What have you learnt?

The blade becomes red hot on heating. On cooling, the blade regains its original colour. This shows that shaving blade undergoes a physical change on




Fig.7.4.2 Heating a shaving blade



heating. Any change in properties such as shape, size, colour and state of a substance is called a physical change. Generally, physical changes are reversible in nature.


Extension

Give a few more instances where change in colour/appearance indicates a physical change.

Notes				

PHYSICAL AND CHEMICAL CHANGES ______12





A chemical change involves formation of a new substance

How will you proceed?

- 1. Take a 5 cm long magnesium ribbon and clean it with a sand paper.
- 2. Hold it with a pair of tongs over a flame and burn.
- 3. What is the colour of flame when magnesium ribbon burns? What is left behind when magnesium ribbon is completely burnt?
- 4. You can see that a white powdery ash is obtained after the burning of the ribbon.





Fig.7.4.3 Burning of magnesium ribbon



- 5. Is it a new substance?
- 6. Collect the white powdery ash and put it in a test tube.
- 7. Mix it with a small amount of water and stir the mixture well.
- 8. Test whether this solution is acidic or basic by adding a drop of phenolphthalein to the solution.
- 9. Do you see any colour change?
- 10. What do you conclude?

Burning of magnesium ribbon forms a new substance Magnesium oxide. This shows that burning of magnesium ribbon is a chemical change.

Extension


Represent the above chemical change in the form of a chemical equation.

PHYSICAL AND CHEMICAL CHANGES _

121







ACTIVITY 7.4.4

To study the nature of change What is required? when zinc is added to copper copper sulphate, sulphate solution.

test tube, zinc granules, spatula.

How will you proceed?

- Dissolve a spatula full of copper sulphate in 2 mL of water in a test tube.
- 2. Add 1 or 2 zinc granules into it and keep it for some time.
- 3. Do you observe any change in colour of the solution and zinc granules?
- 4. You will observe that the zinc granules turn brown and after sometime, the solution turns colourless.
- 5. What may be the possible reason for change in colour?
- 6. The change in colour of zinc granules is due to the deposition of copper on it. The change in colour of the solution is due to the formation of a new substance, zinc sulphate which is colourless.




Fig.7.4.4 Studying the nature of change in a test tube when zinc is added to copper sulphate solution



The reaction between zinc granules and copper sulphate involves formation of two new substances, i.e., zinc sulphate and copper (brown). Thus, this process is a chemical change.


Extension

Write the reaction between zinc granules and copper sulphate in the form of a chemical equation.

Notes			

PHYSICAL AND CHEMICAL CHANGES ______12





1. Take a W-tube and add a pinch of baking soda from

the left side tube.

- 2. Take a few drops of limewater in a dropper and add to the right side of the W-tube.
- 3. Take a few drop of vinegar (acetic acid) in the dropper and fix the dropper at the mouth of the left side tube.
- 4. Why the tube is blocked with the dropper?
- 5. Start adding vinegar drop wise to the baking soda.
- 6. Do you see bubbles of gas coming out?
- 7. What happens to the lime water?
 You will notice that the lime water turns milky.





Fig.7.4.5 Chemical changes in W-tube



What	have	you i	learnt?
------	------	-------	---------

- 1. The reaction between vinegar and baking soda produces carbon dioxide gas which when passed through lime water, turns it milky.
- 2. This reaction involves the formation of new substances. So, it is a chemical change.

Notes			

PHYSICAL AND CHEMICAL CHANGES ______1



Activity 7.4.6

How do you get the crystals Beaker, copper of copper sulphate from its sulphate powder, dil. solution?

How will you proceed?

Take 5 mL of water in a china dish.

G-clamp china dish. 2. Add sufficient amount of copper sulphate powder to get a saturated solution.

3. Add a few drops of dil. H₂SO₄ to this solution.

4. Heat the solution on the kerosene burner and continue adding copper sulphate till no powder more dissolved.

What is required?

H,SO, water, glass

funnel, tripod stand,

rod, filter paper,

kerosene burner,

5. Allow the solution to cool for sometime.






Fig.7.4.6 Arrangement to show the crystalisation



- 6. Filter the contents of china dish in a beaker to get a clear solution.
- 7. Keep the solution for cooling for some time without disturbing it.
- 8. What do you observe after sometime?

Crystals	of copper sulphate are obtained by
•	zation of the copper sulphate solution.
Crystals	of pure substances can be obtained from
their sol	utions. This process is called crystallisation

PHYSICAL AND CHEMICAL CHANGES _______133





7.5 WIND, STORMS AND CYCLONE

Activity 7.5.1

How can we show that air expands on heating and contracts on cooling?

What is required?

A test tube, a rubber cork with one hole, a long fine glass tube, a rubber tube, some coloured water, a wax candle, match box, beaker and water.

How will you proceed?

- 1. Take a test tube.
- 2. Fit a one-holed cork on its mouth.
- 3. Insert tightly a fine glass tube through this hole [Fig. 7.5.1]. Put a little molten wax through the gap in between the cork and glass tube to make the assembly air tight.
- 4. Insert a drop of coloured water ink, KMnO₄ through the glass tube.
- 5. Now hold the test tube in the bright Sun or heat it gently at the end using a candle flame.




Fig.7.5.1 Arrangement to show expansion of air

- 6. Observe the motion of the water drop in the glass tube. What happens to it?
- 7. Now keep the test tube in shade.
- 8. Observe the motion of water drop in the glass tube again.
- 9. What happens to the water drop if the test tube is placed in ice cold water?

Air expands on heating and contracts on cooling.

Extension

Make a slight alteration in the above arrangement. Instead of taking a drop of coloured water in the glass tube, connect one end of the rubber tube to the open end of the glass tube. Immerse the other end of the rubber tubing in water taken in a beaker. Now repeat the activity. Observe the beaker. What happens now? Why?

WIND, STORMS AND CYCLONE .

135







How will you proceed?

- 1. Take a very thin walled plastic bottle.
- 2. Pour small quantity of hot water into it.
- 3. Shake it well.
- 4. Pour out the hot water and close the lid of the bottle immediately.
- 5. Now pour cold water over the bottle.
- 6. Observe what happens.
- 7. Try to find the reason for your observations.
- 8. Now open the lid of the bottle.
- 9. What happens now?



Fig.7.5.2 Bottle with hot water deformed when cooled





How can we show that the air pressure gets reduced when air blows?

What is required?

A funnel, a plastic ball.

How will you proceed?

- 1. Take a funnel.
- 2. Blow through the stem keeping the funnel vertically upwards.
- 3. Now take a plastic ball and place it inside the funnel.
- 4. Now blow air through the funnel.
- 5. Does the ball rise up?
- 6. Now invert the funnel.
- 7. Hold the plastic ball in your hand just below the mouth of the stem of the funnel.
- 8. Blow air through the stem of the funnel and remove your hand.
- 9. Does the ball fall? Try to find the reason?




Fig.7.5.3 Arrangement to show air pressure

WIND, STORMS AND CYCLONE __



What have you learnt?
The air pressure decreases at the mouth of a nozzle
when air is blown though it.
Notes



7.6 Soils

ACTIVITY 7.6.1

How will you show that the soil collected from different places differ in composition?

- Soils are classified on the basis of humus content and particle size.
- The various types of soils are sandy, clayey and loamy soil.

What is required?

Soil samples (from garden, fields), watch glass, glass rod and tumblers(50 ml beaker.

How will you proceed?

- 1. Take two tumblers, each filled 2/3rd with water.
- 2. Add about 10 spoons of soil samples separately in each of them.
- 3. Stir them well with a glass rod.



Fig. 7.6.1 Different layers of soil



- 4. Leave these tumblers undisturbed for sometime. Do you see different layers [Fig. 7.6.1].
- 5. Are the layers similar in the tumblers?


Soils from different places contain different amounts of gravel, sand, clay and humus.

Extension

Repeat the experiment with different soil samples.

Notes			





- 1. Take some soil in a borosil test tube and heat it on a burner with the help of a holder.
- 2. Do you observe water vapours? Do you see water droplets on the inner side of the test tube?




Fig. 7.6.2 Water is present in soil sample

Soils contains water. On heating, the water of the soil gets evaporated and soil becomes loose.

Extension

Carryout the above activity with some more soil samples and compare their water contents.

SOIL ______14



ACTIVITY 7.6.3

How will you show that different types of soils have different water holding capacity?

How will you proceed?

1. Take equal amounts of the three types of soil, viz., sandy, clayey and loamy. What is required?

Plastic funnels, beakers,

filter paper, dropper,

and measuring cylinder.

- 2. Dry them and keep in separate funnels A,B and C lined with wet filter paper.
- 3. Place these funnels on measuring cylinders.
- 4. Gradually pour 100 mL of water in each funnel.

5. After about one hour, read the level of water in the cylinder.




Fig.7.6.3 Water holding capacity of clayey, sandy and loamy soil

6. Record your observations in a table as given below:

	Set up A (with clay)	Set up B (with sand)	Set up C (with loam)
Volume of filtered	mL	mL	mL
water			



Conclusion

- 1. Funnel A has passed......mL water (fill in the blanks by actual reading and calculations) and retained.....mL. The amount of water retained is more than that has passed. Thus, we can say that clay has more water-holding capacity.
- 2. Funnel B has passedmL of water at a faster rate due to the large size of particles and bigger air spaces in between. Thus, the sandy soil has very little water retention capacity (.....mL).
- 3. Funnel C has passed more water than in A and less than in B due to mixture of sand and clay. It also contains humus which retains water. It has retained a moderate quantity of water (.....mL).

What have you learnt?

Texture of soil depends upon the size and arrangement of soil particles. Soils can be classified into different types on the basis of size of soil particles and spaces present in-between which decides its water holding capacity.

Extension

Repeat the above activity with other type of soil and calculate the percentage of water absorbed.

SOIL _______ 143





7.7 RESPIRATION IN ORGANISMS

ACTIVITY 7.7.1

How to demonstrate the presence of CO in the air we breathe out?

The air we breathe in is a mixture of gases. The air we breathe out contains a higher amount of CO₂ (approx.4.4 percent and water vapours) as compared to the air we breathe in (approx.0.04 percent). Respiration is a cellular process resulting in the liberation of energy, whereas breathing is simply the exchange of gases between the organism and the atmosphere.

What is required?
Two test tubes, a
Y-shaped tube, two
rubber corks each
with a hole, lime
water and tap water.

How will you proceed?

- 1. Take two test tubes.
- 2. Take some freshly prepared lime water, in one test tube.
- 3. Take the same quantity of tap water in second test tube
- 4 Insert the two ends of the Y-shaped tube into the test tubes so as to dip the ends the limbs in the liquids.
- 5. Now, blow out gently through the median limb of Y-shaped tube a few times.
- 6. Do you observe any change in the colour of lime water and tap water in the two test tubes?






Fig.7.7.1 Carbon dioxide is present in the air we breathe out

When we blow into the lime water, it turns milky.

Can you give reasons?

Extension

- 1. Blow out a few times on the surface of a mirror. Now, release some air on the surface of another mirror through cycle pump. Did you find any difference? Give reasons.
- 2. Whenever we are sleepy or in a state of low activity, we generally yawn. Can you give the reasons?
- 3. During severe winters, when we exhale, a mist of white smoke comes out from our mouth. Find out the reasons.
- 4. In a crowded and closed area, we often feel suffocated. Why?

RESPIRATION IN ORGANISMS _____

145



ACTIVITY 7.7.2

What is the role of diaphragm in breathing?

What is required?

A glass or plastic 'Y' tube, bell jar with air tight cork having a hole at the centre, two small balloons, a thin rubber sheet or a big balloon.

Breathing means inhalation and exhalation of air. Lungs are the main breathing organs which are situated in the chest cavity surrounded by the ribs. A large muscular sheet, the diaphragm is present below the lungs which forms the floor of the chest cavity. Movement of the diaphragm helps in breathing.

How will you proceed?

- 1. Take a plastic or glass 'Y' tube.
- 2. Fix two deflated balloons to the forked ends of the 'Y' tube.
- 3. The 'Y' tube is kept in inverted position and is inserted into the hole of the cork from top side.
- 4. Fix the cork in the mouth of the bell jar.
- 5. Tie a thin rubber sheet or a split balloon is fited on the open end of the bell jar with the help of a rubber band.
- 6. Pull the rubber sheet down from the base downward, you will see that the balloons at the forked ends swell up.
- 7. Push the rubber sheet up. The balloons will shrink.
- 8. The balloons in the apparatus represent lungs and the bell jar, the chest cavity. The rubber sheet represents the diaphragm.








Fig.7.7.2 Model to show mechanism of breathing

During inhalation, ribs move up and outwards, and diaphragm moves down to increase space resulting in the increase of volume in the chest cavity and air rushes into the lungs. The lungs get filled with air.

During exhalation, the ribs move down and inwards while diaphragm moves up. This reduces the volume of the chest cavity and increases the pressure on the lungs, with the result that the air is pushed out of the lungs.

Extensions

- 1. With the help of a measuring tape, measure the chest size of your classmates during breathing.
- 2. Observe the breathing in fishes and compare it with other organisms.

RESPIRATION IN ORGANISMS _

147



7.8 TRANSPORTATION IN ANIMALS AND PLANTS

ACTIVITY 7.8.1

How can you listen to your own heart beat?

What is required?
Funnel (6-7 cm in diameter), rubber tube (50cm), rubber sheet or balloon.


You can feel your heart beat if you place your one hand on the left side of your chest. But can you hear it? The stethoscope is a device which amplifies the sound. Let us construct our own stethoscope.

How will you proceed?

Rubber tube

1. Take a small funnel and fix the rubber tube tightly on the stem of the funnel.

Stretch a rubber sheet or a balloon on the mouth of the funnel and fix it tightly.



Funnel



3.	Put the open end of the tube on one of your ears and place the mouth of the funnel on your chest over the heart. Do you hear any sound?
W	hat have you learnt?

The regular thumping sound that you hear is the heart beat. Normally the heart beats 72 times a minute in an adult.

Notes			

TRANSPORTATION IN ANIMALS AND PLANTS ______ 149



ACTIVITY7.8.2

How does water move from one cell to another in plants?

What is required?

Potato, 3 per cent sugar/
salt (3gm in 100 mL)

solution, petridish, knife
or scalpel and pin.

How will you proceed?

- 1. Take a large potato and peel off its outer skin.
- 2. Give two cuts on opposite sides to make two flat bases.
- 3. Make a deep and hollow cavity on one of the flat surface leaving only a thin base.
- 4. Fill half of the cavity with 3 per cent sugar solution and mark the level by inserting a pin in the wall of the potato from inside the hollow cavity.
- 5. Put the potato in a petridish containing some water and the level of water should be below the level of the pin.
- 6. Allow the apparatus to stand for a few hours.
- 7. You will observe an increase in the level of sugar solution inside the cavity of the potato.





Fig 7.8.2 Water moves through different cells of potato through osmosis



Water moves from a region of lower sugar concentration to the region of higher sugar concentration across a membrane. By this method plants absorb water and minerals by the roots. The roots have root hair. The root and hair absorb the water present in the soil and transport it to other parts of the plant.

Extension

Take thick peels of potatos, cucumber or any other vegetable. Put them in water. The peels will curve showing absorption of water by cut end.

TRANSPORTATION IN ANIMALS AND PLANTS _____







7.9 Reproduction in Plants

ACTIVITY 7.9.1

How does asexual reproduction occur in plants?

What is required?
Yeast powder, 10 per cent
sugar (10 gms sugar in
100mL of water) solution,
glass beaker, glass rod, slide,
cover slips, microscope.

In asexual reproduction, plants give rise to new plants without seeds. Only a single parent is involved. Budding is one of the process of asexual reproduction.

How will you proceed?

- 1. Fill one-fourth of a beaker with 10 per cent sugar solution and add a pinch of yeast powder or one half powdered yeast tablet to it.
- 2. Keep it in the warm place of the room for an hour.
- Now put a drop of this liquid on a glass slide and put a coverslip. Observe under the microscope under high power.
- 4. You will observe small bulb-like projections coming out of the yeast cells. These are called buds.
- 5. The buds gradually grow and may get detached from the parent cell or result in a chain of buds.
- 6. The new yeast cell grow, mature and produce more yeast cells.





- 1. The process you observed is called budding.
- 2. Budding is a kind of asexual reproduction. In this process single parent is involved and new plants are produced without seeds.




Fig. 7.9.1 Asexual reproduction in yeast by budding


Extension

Observe the methods of propagation of potato, sugarcane, ginger bryophyllum, etc.


REPRODUCTION IN PLANTS _

153





7.10 MOTION AND TIME



How will you proceed?

- 1. Fix the G-clamp at the edge of a table as shown in Fig.7.10.1. Insert the iron rod in the hole of the G-clamp and tighten the screw.
- 2. Tie one end of the thread to the hook of the bob and suspend the bob from the rod of G-clamp making a knot in the thread.




Fig. 7.10.1 Oscillation of a simple pendulum



- 3. Let the bob come to rest in its mean position say at O. Mark this position on the ground as shown in Fig.7.10.1.
- 4. Displace the bob slightly to one side from its mean position and then release it gently.
- 5. Measure the time for 20 oscillations with the help of a stop watch or wrist watch.
- 6. Divide this by 20, to get time for one oscillation.
- 7. Repeat the steps No. 5 and 6 at least three times by displacing the bob slightly more i.e. by increasing 2 cm each time from the mean position.
- 8. Does the time period of the pendulum depend on the displacement of the bob from its mean position?

Time period of a given simple pendulum is always the same and it does not depend on the displacement of the bob from its mean position provided this displacement is not very large.

Extension

Repeat this activity by changing the length of the thread say, from 50 cm to 150 cm and determine the time period in each case.

Is the time period same in both the cases?


If possible, change (a) the diameter, (b) mass of the bob and find the factors on which the time period of a simple pendulum depends.

MOTION AND TIME _____










How will you proceed?

- 1. Mark a straight line of length 50 m in an open space with the help of a measuring tape.
- 2. Ask one of your friend to start walking on this line with his normal gait/walk and simultaneously start the stop watch.
- 3. Stop the stop-watch at the moment when your friend reaches the other end of the straight line.
- 4. Measure the time taken by your friend to travel the distance of 50 m.
- 5. To find the speed of the normal walk of your friend divide 50 m by the time taken(in seconds) to travel this distance.
- 6. Repeat this activity (steps 2 to 5) with as many friends as possible.
- 7. Find the average of the above speeds to determine the average speed of a normal walk.

What have you learnt? Average speed can be determined by dividing the total distance travelled with the total time taken to travel it.







How will you proceed?

- 1. Place the compass needle in an open space where you can get sunshine throughout the day.
- 2. Mark a line along N-S direction.
- 3. Replace the compass needle by the sun-dial.
- 4. Rotate the sun-dial in such a way that its gnomon is along the N-S direction.
- 5. Note the time on the sun-dial scale where the shadow points.
- 6. Compare this time with your wrist watch.

What have you learnt?

Size and position of the shadow of an object can be used to measure time.




Fig. 7.10.1 Oscillation of a simple pendulum

MOTION AND TIME _______ 157



Extension

For exact measurement of time the gnomon of the sun-dial should be kept in geographical N-S direction. Here we have placed the sun-dial in the magnetic N-S direction.


The geographical N-S direction is inclined with the magnetic N-S direction with an angle which varies from place to place.

Try to locate the geographical N-S direction at your place.

Notes			



7.11 ELECTRIC CURRENT AND ITS EFFECT



How will you proceed?

- 1. Connect the resistor between two terminals of the generator using connecting wires having crocodile clips [Fig.7.11.1].
- 2. Rotate the coil of the generator at a moderate speed for about 30 seconds.
- 3. Check the hotness of the resistor just by touching it.




Fig. 7.11.1 Generator connected to a resistor

- 4. Replace the resistor by another resistor and pass current through it for almost the same time by rotating the generator.
- 5. Check the hotness of the resistor again.
- 6. Repeat this activity with two or three other resistors. What conclusion can you draw from this activity.

What have you learnt?	
When electric current passes through a wire , it	
produces heating effect.	

Notes			







ACTIVITY 7.11.2


How can we show that an electric current flowing through a conductor produces magnetic effect around it.

How will you proceed?

- 1. Connect the metallic wire to a dry cell with a switch using connecting wires having crocodile clips at the two ends.
- 2. Place the compass needle near to the metallic wire [Fig.7.11.2].
- 3. Now close the switch and observe the deflection of the compass needle.
- 4. Open the switch and again observe the needle of the compass whether it deflects again to come back to its initial position or not.

What is required?

A metallic wire, dry
cell/generator, switch,
magnetic compass,
connecting wires having
crocodile clips at its
ends and sand paper.



ELECTRIC CURRENT AND ITS EFFECT _

_ 16′



5	Reneat	stens 3	and 4	two	times	and	draw	conclusion.
J.	INCPCat	SICPS 3	anu 4	two	unics	anu	uraw	COLICIUSIOII

6.	The deflection of the needle shows that current
	carrying conductor acts as a magnet.

When electric current passes through a wire, it produces magnetic effect.

Extension

Repeat this activity using generator in place of dry cell.

Notes			





ACTIVITY 7.11.3

How can we make an electromagnet?

What is required?

A long iron nail,
insulated copper wire,
a cell/generator, iron
filings, compass needle,
connecting wires
having crocodile clips
at the two ends,
switch and sand
paper.

How will you proceed?

- Take a long iron nail.
- 2. Check whether it is already a magnet or not, by dipping it in iron filings. Make sure that it is not a magnet.
- 3. Now wind an insulated copper wire over the nail.
- 4. Remove the enamel at the ends of the insulated copper wire using sand paper.
- 5. Connect the ends of the insulated wire to a switch and then to the cell or generator [Fig.7.11.3].
- 6. Now dip the nail in iron filings. What do you observe?
- 7. Do the iron filings cling to the tip of the nail?
- 8. Now switch off the current. Are the iron filings still clinging to the nail in the same way?




Fig.7.11.3 An electromagnet

ELECTRIC CURRENT AND ITS EFFECT



9.	times. What conclusion can be drawn by this activity?
W	hat have you learnt?
	When current is passed through a wire wrapped on

When current is passed through a wire wrapped on an iron nail, the nail becomes a magnet. This arrangement is called electro magnet.

Extension

Repeat the above activity by increasing the number of turns, wound on the nail. Find out how the strength of an electromagnet changes with the number of turns wound on the nail.

Notes			



7.12 LIGHT

Activity 7.12.1


How can we identify different types of lenses and mirrors?

What is required?

Plane mirror, convex mirror, concave mirror, convex lens, concave lens, connecting wires with crocodile clip at both ends and sand paper.

How will you proceed?

- Place the given mirrors and lenses on the table and observe them one by one. On the basis of their look and properties classify them into following two groups:
 - (a) Polished/shining/opaque surfaces having same thickness throughout. These are mirrors.
 - (b) Unpolished transparent, curved surfaces having different thickness. These are lenses.
- 2. Take the mirrors and examine their shining surfaces.
 - (a) If the shining surface is plane, it is a plane mirror.
 - (b) If the shining surface is curved outside, it is a convex mirror.
 - (c) If the shining surface is curved inside, it is a concave mirror.
- 3. Now take the lenses and examine them.
 - (a) If both the surfaces are curved outside, it is a double convex lens.



(b) If both the surfaces are curved inside, it is a double concave lens.




Fig.7.12.1 Different types of mirrors or lenses

What have you learnt?

Mirrors and lenses can be identified by observing them.

Extension

Connect ray streak apparatus to the generator with the help of connecting wires having crocodile clips. Place this arrangement on a white sheet. Rotate the handle of the generator to glow the bulb. When the bulb glows a parallel beam of light is seen in the ray streak apparatus.

Place the mirrors one by one in the path of the parallel beam and observe the nature of the reflected beam in each case.



- If the reflected beam is parallel it is a plane mirror.
- If the parallel beam after reflection from the polished surface meets at a point (converges) in front of the surface, it is a concave mirror.
- If the parallel beam after reflection from the polished surface, appears to diverge from a point behind the surface (mirror), it is a convex mirror.
 - Place the lenses in the path of the parallel beam of light and observe the nature of the beam after passing through the lens.
- If the parallel beam of light after passing through the lens meets (converges) at a point on the other side of the lens, it is a convex lens.
- If the parallel beam of light after passing through the lens appears to come from a point on the same side as the incident beam, it is a concave lens.

I ICUT



Activity 7.12.2

How can we find the nature of the image formed by a concave mirror?

What is required?
Concave mirror,
candle, match box
and thick sheet of
white paper.

How will you proceed?

- 1. Light a candle.
- 2. Fix a concave mirror on the stand and place its reflecting side at a distance of about 25 cm from the candle. Ensure that the mirror surface is perpendicular to the ground.
- 3. Hold a thick sheet of white paper (or a screen) on the other side of the candle (behind the candle) keeping its plane parallel to the mirror surface.
- 4. What do you see on the screen?
- 5. If nothing is visible except a dim patch of light, move the screen a little forward or backward in line with the concave mirror and thecandle, till you obtain on it a clear image of the flame.
- 6. Now shift the candle towards the mirror to a distance of about 2-3 cm and again obtain the image of the candle flame on the white screen as done previously.





Fig7.12.2 Nature of image formed by using a concave mirror



- 7. Repeat step No. 6 three to four times by shifting the candle towards the mirror and observe the image in each case.
- 8. Again repeat the activity by shifting the candle away from the mirror and observe the image on the screen in each case.
- 9. What inferences could you draw from the above observations?
- 10. Now on the basis of your observation try to answer the following question .
 - Whether the image is always obtained on the screen?
 - Whether the image formed is always (a) enlarged and; (b) diminished in size?
 - Whether the image formed is always inverted?

The image due to a concave mirror can be:

- (a) real or virtual
- (b) enlarged or diminished
- (c) erect or inverted

LIGHT ______ 169







Activity 7.12.3

How can we find the nature of the image formed by (a) convex lens (b) concave lens?

How will you proceed?

- Fix a convex lens on the stand and place it on the table in such a manner that its surface is perpendicular to the table.
- 2. Light the candle and keep it in front of the lens.
- What is required? Concave lens, convex lens, candle, match box, lens-stand, thick sheet of white paper (or a screen).
- Hold a sheet of white paper (or a screen) on the other side of the lens. What do you see on it?
- If nothing is visible, move the sheet forward or backward in line with the lens and the candle until you obtain a sharp and clear mage of the flame on the screen.
- 5. Now shift the candle towards the lens to a distance of about 2 cm and again obtain the image of the candle flame on the screen.




Fig7.12.3 (a) Nature of the image formed bye using a concave lens






Fig. 7.12.3(b) Nature of the image formed by using a convex lens


- 6. Repeat step No. 5 three to four times by shifting the candle towards the lenses and observe the image in each case.
- 7. On the basis of your observations, Now try to answer the following questions.
 - —Whether the image is always obtained on the screen?
 - —Whether the image is always erect or inverted?
 - —Whether the image is always enlarged or diminished?
- 8. What inferences do you get from the above observations?
- 9. Repeat the activity with a concave lens and make inferences.

- 1. Nature of the image depends on the distance between the object and the lens in case of a convex lens.
- 2. What ever may be the position of the object with respect to the lens it is not possible to get the image on the screen in the case of a concave lens.

IGHT _______ 171







- 1. Take the white card sheet.
- 2. Draw a circle of 10 cm diameter on it and cut a circular disc out of the card sheet.
- 3. Divide the circular disc into seven equal sectors/parts using the scale.
- 4. Use the water colour to colour the sectors in the order violet, indigo, blue, green, yellow, orange and red.
- Make a hole at the centre of the disc such that a pencil head/axle can pass through it. The disc is now ready for use.
- 6. Sharpen the pencil and insert its tip into the hole such that the pencil protrudes partly out of it and the coloured face of the disc is verticals.





Fig.7.12.4 (a) A disc with seven dours (b) It appears white by rotating

- 7. Hold the pencil from the top with its tip towards you. Now rotate it in daylight so that the disc attached to it also rotates rapidly. What do you observe?
- 8. What is the colour of the disc while it rotates?
- 9. Alternatively, what conclusion can you draw from this activity?

White light is composed of seven colours.				
Notes				

LICHT 17





What is required?
Plane mirror,
transparent tumbler
and water

How will you proceed?

- 1. Take a plane mirror of suitable size.
- 2. Place it in a tumbler as shown in Fig.7.12.5.
- 3. Fill the tumbler with water.
- 4. Place the tumbler near a window such that the direct sun rays falls on the mirror.
- 5. Adjust the position of the tumbler such that the reflected light from the mirror falls on a nearby wall.




Fig.7.12.5 Splitting of the sunlight into its constituent colours



- 6. Place a white sheet of paper at this place of the wall.
- 7. What colours do you see on the white sheet.
- 8. Now disturb the water and observe the white sheet. Are the colours still there? If yes, why?
- 9. As water settles do you again see the colours? If yes, why?

What I	rave	vou l	learnt?
* * * * * * * * * * * * * * * * * * * *	0000	, , , , ,	

- 1. Sunlight can be splitted into different colours.
- 2. The above arrangement works as a water prism.

Notes			

LIGHT _______178





UPPER PRIMARY SCIENCE KIT MANUAL

FOR CLASS VIII





राष्ट्रीय शैक्षिक अनुसंधान और प्रशिक्षण परिषद् NATIONAL COUNCIL OF EDUCATIONAL RESEARCH AND TRAINING





8.1 Crop Production and Management

ACTIVITY 8.1.1

Do plants need manures and fertilisers for better growth?

What is required?

few seedlings of moong/gram or any
other commonly occurring plants, 3 small pots with soil, manure and urea.

The substances which are added to the soil in the form of nutrients for the healthy growth of plants are called manures and fertilisers. Soil supplies mineral nutrients to the crops. These nutrients are essential for the growth of plants.

Manure: It is an organic substance obtained from the decomposition of plants and animals wastes. It increases the water holding capacity of the soil. It makes the soil porous and helps in exchange of gases. It improves soil texture and increases the number of helpful microbes in the soil.

Fertilisers: They are chemical substances which are rich in particular nutrients. Fertilisers are produced in factories. Urea, ammonium sulphate, super phosphate and potash are some well known fertilisers which provide the plants with NPK (Nitrogen, Phosphorous, and Potassium). Soil testing is done to find out the nutrients required for the soil

How will you proceed?


- Select three seedlings of equal size of moong or gram.
- 2. Take three empty pots and mark them A, B, C.





- 3. To pot 'A' add a little amount of soil mixed with a little cow dung manure.
- 4. To pot 'B' put the same amount of soil mixed with few crystals of urea.
- 5. Take the same amount of soil in pot 'C' without adding anything.
- 6. Now plant the seedlings in them and pour the same amount of water in each pot.
- 7. Keep them in a safe open place and water them daily.
- 8. After 7 to 10 days observe their growth.

- 1. Plants in all the pots do not show similar growth.
- 2. Plants in pots A and B show faster growth compared to plant in pot C.
- 3. Compare the growth of plants in pots A and B.
- 4. Manures improve the soil texture and water



fertilisers on plant growth.





holding capacity and make the soil well-aerated which facilitates the healthy growth of plants.

Extension


- 1. Ask the farmers which fertilisers they are using for which crops.
- 2. When do the farmers add manures and fertilisers during agricultural activities and why?

Notes			

CROP PRODUCTION AND MANAGEMENT _______17



8.2 MICROORGANISMS: FRIENDS OR FOES



Microorganisms can survive in a wide range of environmental conditions such as soil, water, hot springs, ice cold waters, saline water, deserts, decomposing bodies etc, and even in volcanic eruptions. They are also present inside the body of animals including humans.

How will you proceed?

- 1. Collect some moist soil from the field in a beaker.
- 2. Add water to it and let the soil particles settle down.
- 3. Take a drop of water from the beaker and place it on a clean glass slide.
- 4. Put a cover slip over it and observe under the microscope. Are you able to see any moving organisms?
- 5. Now, prepare another slide in the same way by placing a drop of pond water on it. You can collect water in a beaker from the pond in your locality.



The slides show tiny organisms. Does it mean that water and soil are always full of microorganisms?

Do they have any significance?

Fig.8.2.1 (a) *Two Algae*




Fig. 8.2.1 (b) Two types of Protozon

MIRCROORGANISMS: FRIENDS OR FOES _______181



Extension

- 1. Pull out any leguminous plant, such as pea and beans, from the field. Carefully look at the roots of these plants. You will find some nodules. Take one nodule and crush it on a glass slide and stain it with a suitable stain. Observe it under the microscope.

 Does it show any organisms? What are these? Do they have any significance?
- 2. Scrap the surface of a leaf or a stone from a pond with the help of a scalpel and place the material thus obtained on a slide. Put a drop of water on it and observe under the microscope. Try to identify the organisms present by comparing them with the figures given in your book.



ACTIVITY 8.2.2

Are microorganisms an integral part of our life?

What is required?

Stale bread, curd,

spoiled food, rotten

fruits, petri dish,

slides, cover slips

and microscope.

How will you proceed?

- Identify the food items or products at your home which are suspected to contain microorganisms. You can take the help of your teacher.
- 2. Collect these items.
- 3. Now, with the help of your teacher, try to find out
 - the type of microorganisms present in each food item.
 - the role played by microorganisms in each of them.
 - whether it is harmful or not?

Tabulate your findings in the table given below (Two examples are already given).

Sl.No.	Food Items	Kind of microorganism	Role and whether beneficial or harmful
1.	Curd	Bacteria (Lactobacillus)	Converts and curdles the milk into curd,(beneficial)
2.	Spoiled food	Bacteria	Poisoning, (harmful)
3.			
4.			
5.			

MIRCROORGANISMS: FRIENDS OR FOES _______183




Microorganisms play an important role in our lives. Some of them are beneficial and some are harmful.

Extension

- 1. Take some curd and divide it into three parts. Mix one part thoroughly to 250 mL of each of the following: warm milk, boiling milk and cold milk. Observe each mixture after 4 to 6 hours. Note down the observation and try to explain the reasons.
- 2. Can you tell what happens if you add a little curd to warm milk without mixing it well?
- 3. Take a slice of bread and cut it into three pieces. Heat one piece till it becomes dry and keeps it on a dry plate. Now place the second piece in an airtight container and third in an open and moist plate. Leave them undisturbed for 3 to 4 days. Observe the difference and find out the reasons.





What is required?

A test tube, water,

sugar, yeast powder,

aluminum foil or cotton.

Fermentation is a process of conversion of sugar into alcohol. For this process, microorganisms are used on a large scale in different industries.

How will you proceed?

- 1. Fill 2/3rd of a test tube with warm water.
- 2. Add a teaspoon of sugar in it and dissolve.
- 3. Now put half a spoon of yeast powder to this sugar solution of the test tube. Smell the solution.
- 4. Thereafter, cover the test tube with alumnium foil or plug its mouth with the cotton.
- 5. Leave the tube in a warm place for four to five hours.




Fig.8.3.1Action of yeast on sugar solution

- 6. Now uncover the tube and smell the solution.
- 7. Did you find any difference in the smell of the solution? Can you relate this smell to that of any other item?

What have you learnt?


After fermentation, the sugar solution smells different.

This smell is due to the action of yeast on sugar solution.

Yeast converts sugar into alcohol and carbon dioxide.

MIRCROORGANISMS: FRIENDS OR FOES _____





Extension


- 1. Take a plastic bottle and put two teaspoons sugar and half a teaspoon of yeast powder in it. Fill half of bottle with warm water.

 Now fix the balloon (deflated) onto the mouth of the bottle. Observe after half an hour. Are you surprised to see the result?
- 2. Take 100 gms flour and add some sugar and yeast powder in it. Make soft dough by adding warm water and cover it. Leave it for a few hours in a warm place. Now observe the swollen dough.

Notes			



8.3 Synthetic Fibres and Plastics



How will you proceed?

- 1. Take an iron stand with a clamp.
- 2. Tie a cotton thread of about 30 cm length so that it hangs freely from the clamp.




Fig.8.3.1 A thread tie with a small polythene bag hanging from the clamp



- 5. Note the weights, or number of marbles required to break the cotton thread and record it in the table given below.
- 6. Repeat steps 1 to 5 with different types of available threads of same thickness and record your observations in the table given below.

Sl.No	Types of Thread	Total weight or number of marbles required to break the thread
1.		
2.		
3.		
4.		

The weights/number of marbles required to break the thread indicate the strength of the fibre Which thread/fibre did you observe to be strongest?

Notes				







ACTIVITY 8.3.2

Does the synthetic fabric absorb the same amount of water as the natural fabric?

What is required?

10 cm × 10 cm piece of cotton cloth,

10 cm × 10 cm piece of nylon cloth or

any other synthetic cloth, measuring

cylinder, 2 beakers, water.

How will you proceed?

- 1. Take two beakers and fill them with equal amounts of water using a measuring cylinder (Else, graduated beakers may be used).
- 2. Soak the cotton cloth and nylon cloth in these beakers separately for 5 minutes each.
- Take the cloth pieces out of beakers and spread them in the Sun. Note the time taken by each cloth for drying.





Fig.8.3.2 Different fibre absorb different density of water

SYNTHETIC FIBRES AND PLASTICS _______189



- 4. Note the volume of water that remains in the two beakers using the measuring cylinder
- 5. Which fabric absorbed more amount of water?
- 6. Which fabric dried faster?

Synthetic fabric absorb lesser quantity of water as compared to cotton fabric and that is why they dry faster than the cotton fabric.

Extension


- 1. Try the above activity using polyester, silk, terylene and acrylic cloth.
- 2. What fabric is used in umbrellas and why?

Notes			





8.4 MATERIALS: METALS AND NON-METALS



When materials are hammered some of them break up into smaller pieces whereas the others spread into the sheets. On this basis, we can classify materials into metals and non-metals. Materials which break up into pieces are non-metals whereas the others which spread into thin sheets are metals.

How will you proceed?

- 1. Take each of the given materials and beat them with a hammer one by one.
- 2. What do you observe?




Fig.8.4.1 Beating a piece of wire by a hammer



3. Write your observations in the table given below.


Sl.No	Object	Observations	Metal / Non- metal
1.	Copper wire	Gets flattened into a strip	Metal
2.	Coal piece		
3.			
4.			
5.			
6.			

What have you learnt?

- 1. Materials like copper and zinc can be spread into sheets when hammered. Thus, they are metals.
- 2. Materials like coal and pencil lead break up into pieces when hammered. Thus, they are non-metals.
- 3. This property of metals by which they can be beaten into thin sheets is known as malleability.







What is required?

Paper strip, plastic strip, wooden stick, circuit board assembly, copper strip, iron nail, a piece of coal, aluminum wire.

How will you proceed?

- 1. Take a circuit board assembly.
- 2. Place one of the above materials, which are to be tested for coductivity under the tap key.
- 3. Push the key so that the circuit gets completed.
- 4. Does the bulb start glowing?
- 5. Test the other materials in a similar manner and record your observations in the given table.

Sl.No	Materials	Good/Poor conductor
1.	Iron nail	
2.		
3.		
4.		
5.		





Fig.8.4.2 Testing good conductors and poor conductors

MATERIALS: METALS AND NON-METALS _________19




What	have you learnt?
1.	Metals allow the electric current to pass through
	them are called conductors while the non-
	metals do not allow the electric current to pass
	through them are called poor conductors.
2.	Materials like copper wire and iron nail are good
	conductors of electricity.
3.	Materials like wooden stick, paper strip and
	plastic strip are poor conductors of electricity.
- - - -	Votes
-	









How will you proceed?

- 1. Take a pinch of rust in a test tube and try to dissolve it in a little amount of water.
- 2. Has the rust mixed completely with water?
- 3. Shake the mixture well for sometime.
- 4. Take few drops of the mixture in two separate watch glasses and test them separately with blue and red litmus papers.
- 5. Do you find any change in the colours of the litmus papers?
- 6. Do you find red litmus changing to blue?




Fig. 8.4.3 To check the acidic or basic nature of rust

MATERIALS: METALS AND NON-METALS ______ 198



Rust mixed with water turns red litmus into blue.

Hence, rust is a basic oxide. Colour of blue litmus

paper does not change.

Extension

Test the nature of dull green coating found on old copper vessels.

Notes			



ACTIVITY 8.4.4

How will you test the nature of the gas evolved on burning of sulphur?

What is required?

Sulphur powder, kerosene burner, deflagrating spoon, double mouth flask and rubber cork fitted with water jet, litmus paper strips.

How will you proceed?

- 1. Take small amount of sulphur powder in the deflagrating spoon.
- 2. Heat it over a kerosene burner.
- 3. What do you observe?
- 4. Name the gas produced.
- 5. Take a double mouth flask with a water jet fixed on its lower mouth.
- 6. Place one small piece each of red and blue litmus papers in it and fix it as shown in Fig.8.4.3.




Fig.8.4.4 Burning of sulphur powder

- 7. Lower the deflagrating spoon with burning sulphur into the flask.
- 8. Do you see water rushing in through the jet?
- 9. What may be the possible reason?
- 10. Did you observe any change in colour in the litmus papers?


What you have learnt?	
Sulphur on burning forms	sulphur dioxide which is
acidic in nature.	

Notes			











How will you proceed?

- 1. Take a small piece of magnesium ribbon.
- 2. Clean its surface with sand paper.
- 3. Put it in a test tube and add few mL of water.
- 4. Put a piece of red litmus paper in it.
- 5. Heat the test tube over a gentle flame for few minutes.




 $\textbf{Fig.8.4.5} \ \textit{Determining the nature of solution}$



6.	Do you observe any change in the colour of litmus paper? You would have seen that the litmus paper turns blue.
W	hat have you learnt?
	1. Magnesium reacts with water to form magnesium
	hydroxide and hydrogen gas is liberated.
	2. Magnesium hydroxide is basic in nature.
	3. Thus, Magnesium reacts with water to form a
	basic solution which turns red litmus to blue.
	Notes








How will you proceed?

- 1. Take a small piece of zinc in a test tube.
- 2. Add about 2 mL of water to it.
- 3. Put a piece of red litmus paper into it.
- 4. Wait for a few minutes.
- 5. Do you see any change in the colour of the red litmus paper? You would observe that there is no colour change.
- 6. Heat the test tube on a gentle flame of the kerosene burner for few minutes.



MATERIALS: METALS AND NON-METALS _______ 20



7. Is there any change in the colour of the red litmus paper? You would observe that the red litmus paper turns blue.	
What have you learnt?	_
Zinc reacts with water on heating to form a basic	_
solution.	_
Extension Repeat this activity using an iron strip.	
Notes	
02 UPPER PRIMARY SCIENCE KIT MANUA	L





ACTIVITY 8.4.7

How do metals and non-metals react with dilute hydrochloric acid?

What is required?

Dilute hydrochloric acid, test tubes, magnesium, aluminium, iron, copper, sulphur powder, charcoal and kerosene burner.

How will you proceed?

- 1. Take small samples of the above materials in different test tubes and label them as A, B, C, D, E and F.
- 2. Add few of drops dil. HCl to each test tube.
- 3. Do you observe reaction taking place in all the test tubes?
- 4. Heat the test tubes in which there is no reaction in cold.
- 5. If a reaction takes place, then bring a burning matchstick near the mouth of the test tube.
- 6. What do you observe?
- 7. Do you observe that the reaction took place in some of the test tubes and the flame of the match stick goes off with pop sound.




Fig.8.4.7 Reaction of metal and non-metal with dilute hydrochloric acid

MATERIALS: METALS AND NON-METALS ______ 203



- 8. Which gas do you think is coming out during the reaction?
- 9. Repeat the above steps and record your observations for each of the materials in the table given below.

Test tube	Metal/Non-metal	Reaction with dilute hydrochloric acid	
		At Room temperature	Warm
A	Magnesium (ribbon)		
В	Aluminium (foil)		
С	Iron (filings)		
D	Copper (wire)		
Е	Charcoal (powder)		
F	Sulphur (powder)		


- 1. Some metals like magnesium, aluminum and iron react with dil.HCl to produce hydrogen gas whereas less reactive metals like copper and non-metals do not react with dil.HCl.
- 2. Metals like magnesium and aluminum react with dil. HCl at room temperature whereas iron reacts on heating.

Extension

Repeat this activity with dil.H $_2SO_4$ with the same samples and record your observations in a table.








What is required?


Solid sodium hydroxide,
test tubes, aluminum
foil, kerosene burner.

How will you proceed?

- 1. Place 3-4 pallets of solid sodium hydroxide in a test tube.
- 2. Add a small amount of water into it to prepare a solution of sodium hydroxide.
- 3. Drop a piece of aluminum foil in this solution,
- 4. Heat the contents on a flame.
- 5. Do you observe any gas evolving from the test tube?
- 6. Bring a burning match stick near the mouth of the test tube. Did it go off with a pop sound?

hydrochloric acid





What have you learnt?

Aluminum reacts with the base (sodium hydroxide)

to form its salt and hydrogen gas is released.

Extension

Repeat this activity with zinc and copper.

MATERIALS: METALS AND NON-METALS _

. 205



ACTIVITY 8.4.9

What happens when metals are immersed in the salt solution of other metals?

A more reactive metal displaces the less reactive metal from its salt solution. No reaction will take place when less reactive metal is immersed in the salt solution of more reactive metal What is required?
Two beakers, copper sulphate, zinc sulphate, ferrous sulphate, zinc, iron, copper strips.

How will you proceed?

- 1. Take copper sulphate solution in two beakers.
- 2. Put an iron strip in one beaker and a zinc strip in the other.
- 3. Do you observe any change in colour of solution?
- 4. Do you see any deposit formed on the surface of metal immersed?
- 5. Repeat the same procedure for
 - (i) copper and iron strips with zinc sulphate solution and
 - (ii) copper and zinc strips with ferrous sulphate solution.
- 6. Record your observations in the table.

Sl.No.	Metal	Salt	Change in Deposit Displacement		Displacement
		solution	colour of	formed	of metal occurs/
			solution	on surface	does not occurs



The more reactive metals like zinc or iron displace the less reactive metal like copper.

Extension

- 1. Arrange the above metals in the decreasing order of their reactivity.
- 2. Which is the most reactive metal?
- 3. Which is the least reactive metal?

Notes			

MATERIALS: METALS AND NON-METALS _



8.5 Combustion and Flame

ACTIVITY 8.5.1

Are all substances around us combustible?

What is required?

Straw, matchstick, wood, paper, iron nail, stone piece, glass, charcoal, a pair of tongs, a glass rod, kerosene burner.

How will you proceed?

- 1. Light the kerosene burner.
- 2. Using a pair of tongs, hold a piece of straw over the flame.
- 3. What happens to the straw?

When a substance reacts with oxygen to give off heat and, light either as a flame or as a glow, this chemical process, is known as combustion.






Fig.8.5.1 Distinguish combustible and non-conbustible substances



- 4. Record your observations in the table given below.
- 5. Repeat the above procedure with other materials and record your observations in the table given below.

Sl.No.	Material	Combustible	Non-combustible
1.	Straw		
2.			
3.			
4.			
5.			
6.			

- 1. Substances like wood, paper, straw etc. burn rapidly in oxygen to give out heat and light and are combustible substances.
- 2. Substances like stone and glass do not produce heat and light when heated in oxygen and are known as non-combustible substances.


COMBUSTION AND FLAME

209










How will you proceed?

- 1. Fix a candle on a table and light it.
- 2. Keep the bell jar over the candle so that it rests on two wooden blocks, and air can enter the bell jar.
- 3. What happens to the flame?
- 4. Remove the wooden blocks and allow the bell jar to rest on table. Now observe the flame of the candle.
- 5. Now cover the open end of the bell jar with a glass plate. What is the effect on the flame?










Fig.8.5.2 Experiment shows that air is essential for burning



- 1. The candle burns freely when the air can enter the bell jar from below. In this way, hot air escapes from the top and fresh air enters from the bottom.
- 2. The flame of the candle flickers and produces smoke when sufficient air does not enter into the flask.
- 3. The flame goes off when the bell jar is covered completely as fresh air is not available.

Extension

- 1. Take a piece of burning wood or charcoal on an iron plate and cover it with a glass tumbler. Observe what happens to the burning wood or charcoal?
- 2. Can you now explain why blankets are wrapped around the person, whose clothes catch fire?

COMBUSTION AND FLAME

. 211







ACTIVITY 8.5.3

Can we heat water in a paper cup over a candle flame without burning it?

What is required?
Two paper cups,
candle, match
box, water.

How will you proceed?

- 1. Hold an empty paper cup over a burning candle flame. What happens to the paper cup? (plastic coated paper cup should not be used).
- 2. Fill the other paper cup with 50 mL water and heat it over the flame. What happens to this paper cup?
- 3. Do you feel that the water in the cup becomes hot?




Fig.8.5.3 Heating water in a paper cup

What have you learnt?

- 1. The water in the paper cup absorbs the heat and the ignition temperature of paper is not reached.
- 2. Substances burn only if they attain their ignition temperature.

Extension

Try to boil water in the paper cup



ACTIVITY 8.5.4

Is the flame hot throughout?

How will you proceed?

- 1. Fix a candle on a stand/ table and light it.
- 2. Hold a copper tube with a pair of tongs and introduce its one end in the innermost part of the flame close to the wick.
- 3. Bring a lighted match stick near the other end of the copper tube. What do you observe?

What is required?

Candle, matchbox,

copper tube of

approximately

5 cm length and

3 mm diameter, glass

slide, pair of tongs,

thin copper wire.




Fig.8.5.4 Detecting the heat in different parts of the flame

- 4. When the candle flame is steady, hold a clean glass slide horizontally into the middle zone of the flame for about 10 seconds and remove it. What do you observe on the slide?
- 5 Now hold a thin copper wire on the outer region of the flame for about 30 seconds. What happens to the wire?

COMBUSTION AND FLAME _

213



- 1. The inner most part of the flame contains unburnt wax vapours which rise into the tube and form flame on the other part of the tube.

 This part of the flame is the coolest part.
- 2. The middle region of the flame contains unburnt carbon particles which get deposited on the glass slide. Partial combustion of wax occurs here.


 This part of the flame has higher temperature than the inner part.
- 3. The outermost region of the flame is the hottest region. Here the temperature is high enough to heat the copper wire strongly and it starts glowing.

Extension

- 1. What is the colour of the non-luminous region of the flame?
- 2. Which zone of the flame is used by goldsmiths for melting of gold and silver?



8.6 CELL: STRUCTURE AND FUNCTION



Organisms may be unicellular (one celled) or multicellular (many celled). These cells vary in shape and size based on their functions.

How will you proceed?

- 1. Take a permanent slide of *Amoeba* from the biology laboratory of your school.
- 2. Study the slide under the microscope. If required, you can adjust the view through the knob present on the microscope.



Fig.8.6.1 Observing the permanent slides of Amoeba and Paramecium



- 3. Now, in a similar way, study the permanent slides of *Paramecium* ands yeast.
- 4. Do they have shape similar to Amoeba? If not, what do they look like?
- 5. Observe all the slides carefully and draw their sketches in your notebook.


Amoeba does not have any definite shape and looks irregular. Whereas, Paramecium looks slipper-shaped and has a definite shape. Yeasts, on the other hand, may show some bulges on their body. Do you know what these are? Recall from the activity performed on yeasts in Class VII.

Extension


- 1. Observe the permanent slides of nerve cell, blood cells and muscle cells. Do they have similar shape or they look different? How do their shapes help them to carry out their functions? Find out.
- 2. Do you know which part of the cell provides it shape?
- 3. Try to obtain the information about the kinds of shapes, other than the above, which are found in cells of different organisms.







- 1. Take a fleshy leaf of onion bulb.
- 2. Peel off the skin from the inner side with the forceps.
- 3. Place the onion skin in watch glass with water.
- 4. Now add a few drops of methylene blue.
- 5. Cut a small piece of the peeled onion skin using a razor blade.
- 6. Now put a drop of glycerin on a slide and place the piece of the stained peeled onion skin on it. Make sure that the peel is not folded.
- 7. Gently place a cover slip over the peel using a needle.
- 8. Clean the edges of the cover slip with a blotting paper.
- 9. Observe the slide under the microscope.



CELL: STRUCTURE AND FUNCTION -

217



Inder the microscope you will find, rectangular shaped structures arranged like a brick wall.

You will see that each structure is outlined by a cell wall and cell membrane and there is a distinct dark structure – the nucleus within the Cytoplasm.

The small structures which look alike are the basic building blocks of the onion. These structures are called cells. All other organisms are also made up of cells. The major difference between a plant cell and an animal cell is the presence of cell wall in the plant cell.

Extension

You can take a fleshy leaf of other plants and take out their peel to observe the cells. To observe animal cells you can make cheek cell preparation. For this scrap your cheek gently with a tooth pick and place it on a slide. Stain the material with methylene blue. Place a cover slip and then observe under microscope.



8.7 REPRODUCTION IN ANIMALS

ACTIVITY 8.7.1

How does asexual reproduction occur in small animals such as Hydra and Amoeba?

How will you proceed?

- 1. Take permanent slide of *Hydra* showing budding and observe under the microscope.
- 2. Did you observe any bulges coming out from the parent body? What are these?
- What is required?

 Permanent slides of

 Hydra showing budding

 and Amoeba showing

 binary fission and

 compound microscope.
- 3. Now observe the permanent slide of *Amoeba* under the microscope.
- 4. Are you able to see the division of parent cell into two parts?
- 5. Draw the diagrams in your notebook and label them.





Fig.8.7.1 (a) Budding in hydra



Hydra and Amoeba reproduce asexually by the process of budding and binary fission respectively.

Extension

Find out other methods of asexual reproduction found in lower animals.



8.8 Force and Pressure

ACTIVITY 8.8.1

How can we apply a force on an object without bringing it into physical contact with another object?

How will you proceed?

A Force between two bar magnets:

1. Fix a bar magnet on a round disc of thick thermocole sheet with the help of candle wax.

2. Put this thermocole disc in a plastic tub filled with water so that it floats on the surface of water.


- 3. Now slowly bring the north-pole of another bar magnet near one end of the floating bar magnet and see what happens.
- 4. Now bring the south-pole of the same bar magnet near the same end of the floating bar magnet and observe what happens now.
- 5. Repeat steps (3) and (4) using the same or different magnets for few more times and draw conclusion based on your observations.
- 6. Can you conclude that one magnet exerts force on the other without having any physical contact?

What is required?

A pair of bar magnets,
thermocole piece,
candle, match box,
comb, two plastic
straws, two pencils
and a plastic tub.

Are these forces attractive or repulsive or both attractive and repulsive?

- B. Force between two charged straws:
 - 1. Now take two plastic straws.
 - 2. Rub a plastic straw with a sheet of paper and put the charged straw over two pencils lying on the table as shown in Fig.8.8.1.
 - 3. Now rub another straw with a piece of paper and bring it near the first straw (Make sure that two straws do not touch each other) and observe what happens.
 - 4. Repeat this activity at least two more times. Draw conclusion based on your observations.



What have you learnt?

Objects experience magnetic or electrostatic forces without being brought into physical contact.

Extension

Observe the objects falling freely towards the ground. Is there any force acting on the falling objects? If yes, what type of force is it? Is it a contact force or a non-contact force?



ACTIVITY 8.8.2

How can we show that force acting on a smaller area exerts larger pressure?

How will you proceed?

- 1. Fill the tray with sand.
- 2. Fix four pointed nails at the four corners of a wooden hard board piece in such a manner that about half length of every nail remains out equally on either side. Ensure that all the pointed ends of the nails remain on the same side of the hard board.
- 3. Keep this arrangement gently in the sand filled tray and put gently a body of about 500 g weight on the hardboard (Fig. 8.8.2).

What is required?

A piece of hard

wooden board (8 cm

× 8 cm), four identical

nails (each 10 cm

long), a body of about

500 g weight, sand

and a plastic tray

(15 cm × 15 cm).






Fig 8.8.2 Arrangement to show the force acting on different area

FORCE AND PRESSURE -

223

- 4. Observe the penetration of the nails of the board in the sand.
- 5. Now reverse the side of the wooden hardboard and again put the same body of 500 g weight on the hard board.
- 6. Observe the penetration of the nails of the board in the sand again.
- 7. In which case do you find more penetration?
- 8. Repeat this activity for few more times and draw conclusion based on your observations.

Same force exerts more pressure when it acts on a smaller area.

Extension

Repeat the above activity with a brick and a bed of sand on which different faces of the brick are kept separately.

In each case also apply impact of a small hammer from a height of about 25 cm. In which case do you find the penetration of brick in sand is more and in which case less?







8.9 FRICTION

ACTIVITY 8.9.1

How can you measure the weight of a given object using a spring balance?

Spring balance is a device used for measuring the force acting on an object. It determines the weight of an object as an attractive force exerted by the earth on the object. What is required?

Spring balance,

G-clamp a small wooden

block or a small stone,

and a piece of thread.

How will you proceed?

- 1. Hold the spring balance vertically and observe its scale carefully.
- 2. If the pointer of spring balance is not at the zero mark, then adjust it at zero.
- 3. Note any two long successive markings on the scale of the spring balance.
- 4. Note down the number of small divisions in between these two long successive markings.
- 5. Divide the value between these two long successive markings by the number of divisions between the two. This is the value of one small division of the scale of the spring balance.
- 6. Now with the help of a thread, suspend the given wooden block from the hook of the spring balance vertically.



- 7. Observe the position of the pointer on the scale of the spring balance and note its reading.
- 8. This is the weight of the given wooden block in terms of mass.




Fig.8.9.1 Measuring the weight of an object by using spring balance

One can use a spring balance to measure the weight of an object.

Extension

A spring balance basically measures the weight of the object attached to it. However, the spring balances usually available in the market are often calibrated in the units of mass (g/kg). Such balances show the ratio of the weight of an object to the acceleration due to gravity at the place of manufacture of the spring balance. The meaning of the acceleration due to gravity will be discussed in detail in Class IX.



ACTIVITY 8.9.2

How can you show that the friction depends upon the nature of the two surfaces in contact?

How will you proceed?

- 1. Attach the wooden block to the hook of the spring balance.
- 2. Spread the sand paper sheet on the table and place the wooden block on it.
- 3. Try to pull the wooden block horizontally on the sand paper sheet first by applying a small force on the spring balance.
- 4. What do you observe?
- 5. Now increase the force gradually on the spring balance till the block starts sliding.

6. Note the position of the pointer of the spring balance, when the block just starts







Fig.8.9.2 Friction between two surfaces



What is required?

sheet (each of size

 $50 \text{ cm} \times 15 \text{ cm}$), a

spring balance, wooden

Sunmica sheet, a sand

paper sheet and a glass



FRICTION ______ 22



- 7. Now repeat the steps No.5 and 6 first with sunmica sheet and then with glass sheet in place of sand paper.
- 8. In which case you require more force to slide the block?
- 9. In which case you require less force to slide the block?
- 10 Draw conclusion based on your observations.

Friction depends upon the nature of the surfaces in contact. For smooth surfaces, friction is less. For rough surfaces, friction is more.

Extension

Repeat this activity by wrapping the wooden block first with a sheet of paper and then with a polythene sheet. Slide the wooden block every time on the surface of the same sunmica sheet and compare the results.







8.10 SOUND

ACTIVITY 8.10.1

How can we distinguish between a pleasant sound and a noise?

What is required?

Four to six glass

tumbler (or bowls) of

same size, a metal rod,

a pencil, water and one

metallic tumbler(bowl).


How will you proceed?

1. Place 4 to 6 glass tumblers (or bowls) keeping separation among them.

- 2. Fill the tumblers (or bowls) with water to different levels.
- 3. Now strike the tumblers (or bowls) one by one with a pencil repeatedly (Fig 8.10.1) and feel the sound produced by each tumbler (or bowl). What do you observe?
- 4. What can you say about the sound produced by each tumbler (or bowl) in Step 3. Is the sound produced pleasant?
- Now strike the metallic tumbler (or bowl) vigoursly using a metal rod.
- 6. Is the sound produced pleasant even now?

What have you learnt?

Some sounds are pleasant and some are not pleasant.



Glasses **Fig.8.10.1**A typical jaltarang producing pleasant sound



ACTIVITY 8.10.2

How can you show that sound can travel through a string?

How will you proceed?

- 1. Make a small hole at the bottom of each cup.
- 2. Tie a string between them by making knots.
- What is required?

 A string (about 10 m long), two paper cups or plastic cups or empty ice cream cups.
- 3. Stretch the string between the two cups.
- 4. Use one cup as receiver and other as mouth piece.
- 5. Now keeping the string sufficiently stretched ask your friend to speak some thing from the mouth piece.
- 6. Try to hear the voice of your friend through the receiver. Can you hear it clearly [Fig. 8.10.2]?
- 7. If you are able to hear the sound of your friend through the receiver then sound has travelled from the mouth of your friend to your ear through the string only.
- 8. Try to draw conclusions from this activity.






Fig 8.10.2 Sound travels through the string

What have you learnt?

Sound can travel through a string.







What is required?

A slinky.

How will you proceed?

- 1. Hold one end of a slinky lying on a floor.
- 2. Ask your friend to sit at a distance holding the other end of the slinky.
- 3. Let him hold the slinky steadily. Now give lateral (side way) jerks to the slinky continuously and observe the pattern of the disturbances (waves) produced in it.
- 4. Increase the frequency of the jerks and observe the pattern.
- 5. Now give larger jerks to the slinky and observe the disturbance that occurs in the slinky from its original position.
- 6. When do you see larger disturbance in the slinky?
- 7. When do you see smaller disturbance?
- 8. Now let the slinky come to rest and ask your friend to hold it.
- 9. Give the slinky a jerk along its length and observe the pattern of the disturbance produced in it.
- 10. Now give this kind of jerks continuously to the slinky and observe the pattern of the wave produced in it.
- 11. Draw conclusion based on this activity.




Fig 8.10.3 *A slinky*

SOUND _____2



Different types of waves can be produced in a slinky.

Extension

Repeat the above activity using a string.

Can we produce similar type of waves in the string as observed above in the case of slinky?

Notes			



8.11 CHEMICAL EFFECTS OF ELECTRIC CURRENT

ACTIVITY 8.11.1

How can you show that water can be splitted into its constituent elements?

Passage of electric current through many substances in molten state causes their decomposition. This process is called electrolysis. Electrolysis of water can be carried out in laboratory where water decomposes into its constituent elements, namely; hydrogen and oxygen.

What is required?

Bell jar, rubber cork

with two stainless

steel electrodes, water,

6 V battery (four dry

cells of 1.5 V each in

a cell holder), plug

key connecting wires,

glass rod, universal

indicator, sodium

sulphate, dropper.

How will you proceed?

- 1. Take water in a 50 mL beaker.
- 2. Add sodium sulphate to water and stir well to prepare its concentrated solution.
- 3. Add universal indicator drop by drop till colour of the solution becomes dark green.
- 4. Take two plastic graduated droppers of 3 mL capacity and cut their lower ends. Fill the remaining 3 cm of its length stem completely with sodium sulphate solution.



- 5. Fix the rubber cork with two stainless steel electrodes in the bell jar and clamp it in the inverted position as shown in Fig. 8.11.1.
- 6. Fill the bell jar with the remaining coloured solution.
- 7. Lower the droppers filled with coloured solution over

two steel electrodes carefully.

- 8. Connect a 6 V battery to the 2 electrodes with connecting wires.
- 9. Also insert a plug key in the circuit and keep the key open.
- 10. Now close the key and observe the two electrodes carefully. What do you notice?



Fig 8.11.1 Water splitted into its components

- 11. Do you observe any bubbles forming on the two electrodes?
- 12. Do you notice any colour change in the two droppers?
- 13. Find out the reason of the colour change.

What have you learnt?

When an electric current passes through it, water splits into two gases-hydrogen and oxygen.

Extension

Repeat the activity with distilled water without adding any salt and note down your observations.





ACTIVITY 8.11.2

How can you use electric current for coating a metal over the surface of another metal?

How will you proceed?

- 1. Prepare fresh copper sulphate solution by dissolving one teaspoonful of copper sulphate in water taken in a 50 mL beaker.
- 2. Add about 10-12 drops of dil. sulphuric acid to this solution.
- 3. Clean the copper plates with sand paper and rinse them with water.

What is required?

6 V battery (four 1.5 V dry cells in a cell holder), 50 mL beaker, copper sulphate crystals, dil. sulphuric acid, connecting wires having crocodile clips at the two ends, 2 copper plates of size about 5 cm × 2 cm, sand paper and plug key.

4. Connect the two copper plates to the +ve and -ve terminals of a battery using separate connecting wires fitted with crocodile clips .





Fig.8.11.2 A simple electric circuit showing electroplating

CHEMICAL EFFECT ON ELECTRIC CURRENT ______ 23



- 5. Now immerse the two electrodes (i.e. copper plates) in copper sulphate solution keeping them at a distance. Take care that the 2 plates do not touch each other [Fig.8.11.2].
- 6. Allow the passage of current through the solution by closing the plug key. Wait for about 20 minutes.
- 7. Now open the key and remove the plates from the solution and observe them carefully. Do the plates look alike?
- 8. Do you find some coating on any of the plates? Note the colour of the coating and also the terminal to which it is connected.


When an electric current is passed through copper sulphate solution the copper sulphate dissociates and copper from the solution goes to the plate connected to the negative terminal of the battery (cathode). To compensate for this loss of copper the same amount of copper from the anode goes to the solution. Thus copper from the plate connected to the positive terminal gets deposited on the plate connected to the negative terminal. The process of depositing a layer of desired metal on another material using electricity is called electroplating.

Extension

- 1. Chrome plated handles of bicycles are very popular these days. How are they made?
- 2. At times ornaments which appear to be made of gold are available at cheaper rates in some shops. How it happens?







8.12 Some Natural Phenomena

ACTIVITY 8.12.1

How is an electroscope
made?

Bell jar, J-shaped metal
wire passing through
rubber cork and having
metallic disc fitted at

How will you proceed?


- Cut a rectangular piece of aluminium foil of size about 10 mm × 10 mm and fold it in V-shape.
 - of scissors, plastic strip
 in with wollen cloth.

 aluminium foil at the lower end

the other end, comb,

aluminium foil, a pair

- 2. Hang the V-shaped aluminium foil at the lower end of the J-shaped wire fitted into the bell jar as shown in Fig.8.12.1. Your electroscope is now ready.
- 3. Rub a comb with dry hair or plastic strip with wollen cloth.
- 4. First bring the rubbed comb near to the disc of the electroscope and then touch it.



J-shaped wire

Fig.8.12.1 Spreading out of the aluminium wire


- 5. Observe the shape of the aluminium foil. Do you notice any change? The tips of the two leaves of the aluminum foil spread out.
- 6. Find the reason of this spread.

hat have you learnt?	
Similar charges are acquired by the two leaves of	the
aluminium foil and they spread out.	
Notes	










- 1. Inflate two balloons.
- 2. Hang them in such a way that they are at the same height but do not touch each other.
- 3. Rub both the balloons with a woollen cloth and then release.
- 4. What do you observe?
- 5. Why do they move away from each other?
- 6. Repeat the above activity:
 - (i) by taking two used ball pen refills in place of two balloons and polyester/woollen cloth. Place one used refill after rubbing with polyester/woollen cloth in a beaker and bring the other one near to it after rubbing.



SOME NATURAL PHENOMENA _



(ii)	by taking one balloon and one refill. What
	difference do you notice in the above two cases
	and why? You may observe the following:
	In activity (i) two refils will be seen to repel each
	other.
	In activity (ii) balloon and refils will be seen to

7. Draw conclusions based on the above activities.

What have you learnt?

attract each other.

Like charges repel and unlike charges attract each other.

Extension

Perform the above activity using two small thermocol balls after suspending them from a support.

Notes			



8.13 LIGHT

ACTIVITY 8.13.1

How can you verify laws of reflection of light using a plane mirror?

What is required?

Generator, ray streak apparatus, plane mirror, connecting wires with crocodile clips, white sheet of paper and sharp pencil.

How will you proceed?

- 1. Close all the slits of the ray streak apparatus except one of it and place the apparatus on a white sheet of paper.
- 2. Connect the bulb of ray streak apparatus with the terminals of the generator.
- 3. Rotate the handle of the generator and observe the ray of light obtained on the sheet of paper.
- 4. Now place a plane mirror in the path of this ray of light.
- 5. What happens to the direction of ray of light after reflection from the mirror?
- 6. Using a sharp pencil draw the positions of the plane mirror, the incident ray and the reflected ray.
- 7. Repeat this activity for the rays incident at different angles on the plane mirror by rotating the mirror.
- 8. Draw normals at the point of incidence in each case.
- 9. Measure the angle of incidence and the angle of reflection in each case and note down their values in the tabular form.



- 10. What relationship do you notice between the angle
- 11. Now take a sheet of stiff paper or a thin sheet of white paper. Let the sheet of stiff paper project a little beyond the table.
- 12. Cut the projected portion of the paper sheet at the middle.
- 13. Repeat the activity given above.
- 14. Make adjustment in such a way that the reflected ray lies on the projected portion of the sheet of paper.
- 15. Now bend the projected portion of paper, on which the reflected ray falls downwards and observe what happens.
- 16. Do you still get reflected ray on the bent portion of the paper? If not, what conclusion can you draw from this observation?




Fig.8.13.1 Reflection of light on a plane mirror

- Angle of incidence is equal to the angle of reflection.
- 2. Incident ray, reflected ray and the normal to the reflecting surface at the point of incidence, all lie in the same plane.





ACTIVITY 8.13.2

How can you show that a reflected light can be further reflected?

What is required?
Two plane mirrors of size about 15 cm X 20 cm and a chair.

How will you proceed?

- 1. Sit on a chair with a window behind you.
- 2. Are you able to see the objects outside the window? Certainly not.
- 3. Now take two plane mirrors.
- 4. Hold one of the mirrors in your right hand and the other in your left hand.
- 5. Now adjust the two mirrors in such a way that you can see the objects outside the window in the mirror, as shown in the Fig.8.13.2.
- 6. Try to find, "How do the rays from the objects behind you travel to your eyes?"
- 7. If possible, try to trace the path of light from the objects to your eyes.




Fig.8.13.2 Reflection of light


What have you learnt?

Light can undergo multiple reflection.

LIGHT ______ 243







How will you proceed?

- 1. Take three rectangular plane mirror strips and join them in a triangular shape with cellotape such that their polished surfaces are on the inside and inclined at an angle of 60° with respect to each other as shown in Fig. 8.13.3.
- 2. Slip the rubber bands on the triangular tube, so formed, at different positions such that the mirror strips remain in place. Fix this triangular tube in a cylindrical tube made of card sheet. Ensure that this tube is slightly longer than the triangular tube.



Fig.8.13.3 Arrangement for making a kaleidoscope

- 3. Take two triangular glass plates of such dimensions as to fit at one end of the tube. Cover one glass plate with tracing paper .
- 4. Place one of the triangular glass plate at one end of the tube such that it is a little inside the tube.
- 5. Now place some multi-coloured bangle pieces over it.
- 6. Cover this end of the tube with another triangular glass plate covered with tracing paper.
- 7. Cover the other end of the tube with a piece of card sheet and make a hole in the centre so that the inside view can be seen through this hole. The kaleidoscope is ready.
- 8. Peep through the hole directing the other end of the tube towards light.
- 9. What do you observe?
- 10. Now rotate the tube a little and observe. Do you get the same pattern now? Why do we see different coloured patterns every time?
- 11. What do you infer from these observations?

We see different coloured patterns in a kaleidoscope
due to multiple reflections of light.

LIGHT ______24









How can you show that the sunlight can be splitted into different seven colours?

What is required?

A triangular glass

prism and a

white screen.

How will you proceed?

- 1. Hold the glass prism such that one of its faces is towards the sunlight while the other is towards a white screen (white wall or white sheet of paper).
- 2. Adjust the prism such that violet, indigo, blue, green, yellow, orange and red light is seen on the white screen.
- 3. What does this indicate?



Fig.8.13.4 Sunlight splitting into seven constituent colours

What have you learnt?

White light is splitted into seven constituent colours when passed through a glass prism.

